



2016

Product catalogue

SUMMARY

Introduction

6



Part 1 Surgery

14



Part 2 Healing

74




A complete range	6
Prosthetic solutions to suit all cases	8
Treatment longevity	10
Primary stability and reliability	11
Quality guarantee	12
Clinical studies	13

Tapered Implants	16
Naturactis	18
Natural+	20
Instruments & common surgical kit	22
Cylindrical Implants	28
Natea+	30
Aesthetica+ ²	32
Instruments & common surgical kit	34
Uneva+	42
Instruments & surgical kit	44
Narrow Implants	60
Naturactis Ø 3	62
Natural+ Ø 3	64
Instruments & common surgical kit	66
Obi Ø 2.7	68
Instruments & surgical kit	70
Peripheral products	64
Common Instruments	66
Extraction kit	70
Bone filling Macrobone®	72

Healing abutments	76
Iphysio® Profile Designer	80

**Part 3
etk prosthesis**

82



**Part 4
Other prosthesis brands**

178



**Part 5
Services**

186



Standard prosthetic components 82

Naturactis / Naturall+
Natea+ 84
Aesthetica+² 110
Uneva+ 126
Naturactis/Naturall+Ø3 144
Obi 156

Customised prosthesis 160

Titanium Esthetibase interfaces 162
Scanbody 164
CAD-CAM works 166

All in bar[®] system 168

Common instruments 172

Prosthetic torque wrench 174
Prosthetic kit 175
Keys & mandrels 176

Titanium Esthetibase interfaces 180
Scanbody 182
CAD-CAM works 184

Welcome pack 188

Ordering and organisation 188

Organisational documents 188
Documents to help with your order 188

Patient information 189

Patient website 189
Waiting room video 189
Patient card 189
Patient file 189
Patient leaflet 189
Waiting room poster 189

YOUR DENTAL IMPLANT PRACTICE MADE EASY!

Implants and services: unrivalled patented implant solutions to streamline your protocols and treatment plans.

SURGERY

Implants
and instruments

Standard instruments and surgical kits

PROFILE DESIGNER Healing and impression

iphysio®

Patented / 3-in-1 / Compact protocol

SYNTHETIC BONE SUBSTITUTE

Macrobone®

PROSTHESIS

1 prosthetic range
for all implants

CAD-CAM
customised prosthesis

All in bar system*

** Unrivalled, patented loading solution that is truly immediate*

TRAINING theory and practice

Dental surgeons

Dental technicians

Assistants

Training courses and integration solutions

SERVICES

Welcome pack


Patient communication material

Ordering and
organisation tools

A COMPLETE RANGE TO SUIT ALL CASES

A complete range


INTRODUCTION



TAPERED IMPLANTS

- Post-extraction sites: tapered implants suitably fill the space left behind by extracted teeth.
- Suitable for sites with very low bone quality.


p 16



CYLINDRICAL IMPLANTS

- Suitable for sites with bone quality with good density.
- Anterior or posterior mandible areas.

p 28




NARROW IMPLANTS

- Gaps where standard implants are not suitable to obtain the desired aesthetic effect.
- Limited vestibular-palatine/lingual ridge size.
- Insufficient mesiodistal space to fit an implant whilst respecting the intervals to be kept with adjacent elements.
- In some cases, using narrow implants can avoid resorting to bone grafting.

Obi implant to be used only for multi-unit restorations.



SHORT IMPLANTS 6 MM



- Effective alternative in clinical situations where the height of the bone is reduced and therefore not allowing for conventional implantology due to the proximity of at-risk anatomical elements.
- Can simplify protocols and avoid resorting to bone grafting.
- Reliable and more easily acceptable therapeutic solution for the patient: reduced cost and treatment duration, improved post-operative outcomes etc.




PROSTHETIC SOLUTIONS FOR ALL CASES

STANDARD PROSTHETIC COMPONENTS

etk has developed a complete range of pre-manufactured abutments and secondary components, enabling you to treat all clinical cases.

This wide range of parts for cemented, screwed and removable prosthesis enables you to carry out single, multiple or total restorations.

- **Nitrided abutments**
 - Preservation of the biocompatibility of the titanium and gingival integration.
 - Yellow colour: less visible under the ceramics.
- **Laser marking of abutments**
 - Better identification of the emergence profile and supra-implant height.
- **Anti-unscrewing design of prosthetic screws**
 - All screws are treated with anti-unscrewing technology, improving mechanical strength and sealing of the implant/abutment junction.
- **Captive screws**
 - Screws are secured thanks to interior threading of prosthetic parts to avoid unwanted collapse.




CAD-CAM CUSTOMISED PROSTHESIS

teknikalab, our expert production centre designs and manufactures your CAD-CAM prosthesis.

- **100% made in France**
 - Raw materials of certified European origin.
 - Perfect traceability from raw material to dispatched product.
- **Offer to undertake comprehensive work**
 - Wide range of materials used: zirconia, IPS e.max®, titanium, chrome-cobalt and PMMA.
 - All CAD-CAM works on the leading implant brands and natural teeth.
- **Fast turnaround**
 - 24-48 hours for simple works (sleeves and bridge frameworks).
 - 48-96 hours for more complex works (customised abutments, implant bars).
- **Your close partner**
 - Phone support by dental technicians.
 - Group or individual training to develop your CAD-CAM skills.
 - Online ordering available 24/7.

p 160

p 178




TREATMENT LONGEVITY

STAE® surface treatment, developed by etk, is the result of 20 years' clinical experience.

This patented surface treatment guarantees perfect osseointegration of implants and therefore longevity of your treatment.


TITANIUM OXIDE MICRO-SANDBLASTING AND ETCHING WITH NITRIC AND HYDROFLUORIC ACIDS

- Subtractive, non-polluting treatment.
- Bone-implant contact surface increased by 79.2%.
- Increased wettability.
- Stimulates osteoconduction.




SPECIFIC AND MASTERED SANDBLASTING PROCESS

- Constant sanding pressure and speed.
- Constant grit size.
- Creation of a macrostructure.



ETCHING PROCESS/ACID PASSIVATION AND NEUTRALISATION

- Decontamination of the surface and formation of a TiO₂ (titanium dioxide) layer.
- Creation of a microstructure.



Numerous histological studies

• **Much research has been carried out on the surface state** of the implants in collaboration with Professor Chappard from the Histology Laboratory at Angers University (France), Prof. P. Bravetti from the faculty of dentistry in Nancy (France), Dr. Giner from the International University of Catalonia (Spain), Prof. Jabbour from St Joseph Faculty of Medicine in Beirut (Lebanon) and the Faculty of Medical Sciences at the University of Iran.

• **Histological observations** on our implants show a large percentage of bone-to-implant contact between the titanium and the cortical bone, as well as newly formed bone of normal texture.

• **In a comparative statistical study** carried out by Professor Chappard on Brånemark® implants, it was shown that etk implants are as reliable as mainstream implant brands.



• **The study on the osseointegration** of 15 etk implants conducted by the University of Iran and published in the Dental Research Journal (Vol. 8, No. 3, January 2014) shows a B.I.C. of 76.82% and an ISQ of 70.83 after 4 months.

PRIMARY STABILITY AND RELIABILITY



MICROTHREAD SYNCHRONOUS WITH THE MACROTHREAD

- Insertion without tearing the cortical bone.
- Stabilises the cortical bone.
- Optimised primary stability.



ASYMETRICAL THREADING

- Homogeneous distribution of masticatory forces.
- Excellent primary stability immediately after implant fitting.

CENTRAL PROTRUSION BETWEEN THREADS

- Contact surface with bone tissue increased by 15%.
- Promotes osteogenesis.
- Activates cellular reconstruction.

DOUBLE HELIX

- Decreased heating of the bone and insertion time.



ENGAGING AND ANATOMICAL APEX

- Threading beginning at the apex, giving the implant a high self-tapping capacity.
- Safe to use in the sinus floor region.



NATURACTIS: DEEP BLADE DESIGN APEX

- Threads blade design at the extremity.
- Allows for a more suitable choice of axis.


YOUR QUALITY GUARANTEE

Thanks to a 100% European integrated design and production process, **etk** ensures total control of processes, the materials used and production conditions (respect of asepsis and the environment).



- **LIFETIME GUARANTEE FOR IMPLANTS***
- **10 YEAR GUARANTEE FOR PROSTHETIC COMPONENTS***
- **teknikalab SECONDARY COMPONENTS: 5-10 YEAR GUARANTEE***
 - Zirconia: 5 years
 - Chrome-cobalt and titanium: 10 years

 More information



** The guarantee only applies subject to the exclusive use of the etk components during all stages of treatment (surgery, healing, impression and prosthesis) and only if all application conditions are met.*

CLINICAL STUDIES

In order to ensure the reliability of our developments and to evaluate our implant systems, we have always been involved in studies with a variety of different partners from universities in France and around the world.

Clinical study of Naturactis dental implants post-extraction dental procedures

J. Ripollés de Ramón, R. Gómez Font, C. Bascones-Ilundain, J. Bascones-Ilundain, A. Bascones-Martínez - University of Madrid (Spain)

Geriatric narrow implants for wearers of full dentures: clinical aspects and prospects with OBI mini-implants

Cédric Huard, Marion Bessadet, Emmanuel Nicolas, Jean-Luc Veyrone - University of Auvergne (Clermont-Ferrand - France)

Placement of implants in the mandible reconstructed with free vascularized fibula flap: comparison of 2 cases with Aesthetica+ implants

Mehmet Kürkcü, DDS, MSc, PhD, Mehmet Emre Benliday, DDS, Cem Kurtoglu, DDS, PhD, and Erol Kesiktas, MD, Adana - Cukurova University (Turkey)

Contribution of a hybrid synthetic and innovating product in bone surgery and its filling: Matri™ BONE with Natea and Naturall implants

Augusto André Baptista, Pierre Bravetti - Henri Poincaré University (Nancy - France)

Multicentre study on the evolution of 3000 Euroteknika and Nobel Biocare® implants from 1984 to 1997 - Comparison of results

Daniel Chappard - LHEA - Angers Faculty of Medicine (France)

Histology and histomorphometrical comparative study of the Universal implant by Euroteknika

Laboratoires Karl Donath, Hamburg (Germany) - Guy Huré, Laboratoire d'Histologie d'Angers (France)

Implant-supported prosthetic solution in cases of small inter-alveolar distance on Aesthetica+ implants

Victor Degasyuk, Ljudmidia Degasyuk - Polyclinic Kiev (Ukraine)

Quantitative study on the roughness of the surface of titanium dental implants and their microstructures

Bally, Dehmas, Rapin - Henri Poincaré University (Nancy - France)

SR Phonares and OBI mini-implants: a perfect fit - Prosthetic strategy

Y. Gastard (Dental Prothesist), F. Truchot, X. Ravalex, G. Bader

Analysis of the purity of surface treatments of Euroteknika implants and its competitors

Jordi Ferre, Joseph Miquel & Giner - Spanish National Research Council (CSIC) - University of Barcelona (Spain)

A comparison of two types of decalcified freeze-dried bone allograft in treatment of dehiscence defects around implants in dogs

Ahmad Moghareh Abed, Rasool Heidari Pestekan, Jaber Yaghini, Seyed Mohammad Razavi, Mohammad Tavakoli, Mohammad Amjadi - University of Iran

Evaluation of the sealing of the connections of Euroteknika implants

Josep Cabratosa Termes, Zaira Martínez Vargas - University of Catalonia (Spain)

An in vitro study to compare the insertion torque and the removal torque of two screw type dental implants with different thread designs on three different materials

Josep Cabratosa Termes, Zaira Martínez Vargas - University of Catalonia (Spain)

Resonance frequency analysis, insertion torque, and bone to implant contact of 4 implant surfaces: comparison and correlation study in sheep

Maroun Dagher, DDS, CAGS, MScD,* Nadim Mokbel, DDS, MSc, PhD, Gabriel Jabbour, DDS, and Nada Naaman, DDS, PhD (Lebanon)

Study of the sealing between the implant and different abutments


H. Ghandi, P.K. Kimani, I. Abou-Rabii, and E. Lynch, University of Warwick, Coventry (UK)

Comparison of leakage at the implant to abutment connection in several implants types using a gas flow method

M.-A. Fauroux, C. Anxionnata, C. Biensa, M. Mechalia, O. Romieua, J.-H. Torresa, Dentistry Department - CHRU de Montpellier (France)

 Download the technical and scientific





Part 1

Surgery

16 Tapered implants

- 18 Naturactis
- 20 Naturall+
- 22 Instruments & common surgical kit

28 Cylindrical implants


- 30 Natea+
- 32 Aesthetica+²
- 34 Instruments & common surgical kit
- 42 Uneva+
- 44 Instruments & surgical kit

50 Narrow implants

- 52 Naturactis Ø 3
- 54 Naturall+ Ø 3
- 56 Instruments & common surgical kit
- 60 Obi Ø 2.7
- 62 Instruments & surgical kit

64 Peripheral products

- 66 Common instruments
- 70 Extraction kits
- 72 Bone filling Macrobone®




1

Surgery

Tapered implant systems




- 18 Naturactis
- 20 Naturall+
- 22 Instruments & common surgical kit





Naturactis

NATURALLY ACTIVE IMPLANT



-
- Internal hexagonal conical connection
 - the same for all implant diameters
 - compatible with Naturall+ and Natea+ implants

-
- Cylindroconical implant
 - Subcrestal placement
 - High primary stability
 - Strong apical anchoring


INDICATIONS

- Post-extraction surgery
- Immediate loading
- Areas with low bone density (D3-D4)

CONTRAINDICATIONS

- Sinus floor region

Technical characteristics



A	\varnothing 3.4 - 3.8 - 4.3 - 4.8 mm
B	Microthread 3 mm
C	Thread 0.4 mm
D	\varnothing 3.5 - 4 - 4.5 - 5 mm
E	Real screw thread 2.4 mm
F	Thread 1.2 mm
G	\varnothing 1.6 - 1.8 - 2.2 - 2.3 mm
H	Sandblasted and etched length

Implant references

Medical grade titanium.

Implants are delivered with a cover screw.

Length H	\varnothing 3.5 mm	\varnothing 4 mm	\varnothing 4.5 mm	\varnothing 5 mm
6 mm	-	-	NIP 45 060	NIP 50 060
8 mm	NIP 35 080	NIP 40 080	NIP 45 080	NIP 50 080
10 mm	NIP 35 100	NIP 40 100	NIP 45 100	NIP 50 100
12 mm	NIP 35 120	NIP 40 120	NIP 45 120	NIP 50 120
14 mm	NIP 35 140	NIP 40 140	NIP 45 140	NIP 50 140
16 mm	NIP 35 160	NIP 40 160	NIP 45 160	-
18 mm	NIP 35 180	NIP 40 180	-	-



Naturall+

ANATOMICAL IMPLANT




-
- Internal hexagonal conical connection
 - the same for all implant diameters
 - compatible with Naturactis+ and Natea+ implants

-
- Tapered implant
 - Bone level
 - High primary stability

INDICATIONS

- All areas
- All bone densities
- Sinus floor region
- Post-extraction surgery
- Immediate loading

Technical characteristics



A	\varnothing 3.7 - 4.2 - 4.7 - 5.2 mm
B	Supracrestal smooth neck 0.4 mm
C	Thread 0.3 mm
D	Microthread 3 mm
E	\varnothing 3.5 - 4 - 4.5 - 5 mm
F	Real screw thread 1.8 mm
G	Thread 0.9 mm
H	\varnothing 2.6 - 3 - 3.4 - 3.8 mm
I	Sandblasted and etched length

Implant references

Medical grade titanium.


Implants are delivered with a cover screw.

Length I	\varnothing 3.5 mm	\varnothing 4 mm	\varnothing 4.5 mm	\varnothing 5 mm
6 mm	-	-	NICP 45 060	NICP 50 060
8 mm	NICP 35 080	NICP 40 080	NICP 45 080	NICP 50 080
10 mm	NICP 35 100	NICP 40 100	NICP 45 100	NICP 50 100
12 mm	NICP 35 120	NICP 40 120	NICP 45 120	NICP 50 120
14 mm	NICP 35 140	NICP 40 140	NICP 45 140	NICP 50 140

INSTRUMENTS

Drills


Point drill Ø 1.5 - 2.2



The upper part of the point drill has the same diameter 2.2 mm of the following drill to prepare its insertion.
The 90° cutting edge makes a clear and precise mark on the bone crest, even when the crest is thin.
Included in the kit.

Reference CFP 15 22 50


Point drill Ø 2.2



The upper part of the point drill has the same diameter 2.2 mm of the following drill to prepare its insertion.
Included in the kit.

Reference NFP 22 180

Initial drills Ø 2.2




Included in the kit.

Lengths	References
6	AFI 22 060
8	AFI 22 080
10	AFI 22 100
12	AFI 22 120
14	AFI 22 140
16	AFI 22 160
18	AFI 22 180

Drills - continuation


Staged drills



Included in the kit.

Diameters	Lengths	References
Ø 2.2 - 2.8	14	NFP 22 28 140
	18	NFP 22 28 180
Ø 2.8 - 3.3	14	NFP 28 33 140
	18	NFP 28 33 180
Ø 3.3 - 3.8	14	NFP 33 38 140
	18	NFP 33 38 180
Ø 3.8 - 4.3	14	NFP 38 43 140
	18	NFP 38 43 180
Ø 4.3 - 4.8	14	NFP 43 48 140
	18	NFP 43 48 180


Naturall+ cortical drills



Included in the kit.

Diameters	References
3.5	NFE 35 22
4	NFE 40 26
4.5	NFE 45 30
5	NFE 50 34


Naturall+ very hard bone drills



Included in the kit.

Diameters	References
3.5	NFC 35 DO
4	NFC 40 DO
4.5	NFC 45 DO
5	NFC 50 DO


Naturactis taps



Included in the kit.

Diameters	References
3.5	NTP 35 100
4	NTP 40 100
4.5	NTP 45 100
5	NTP 50 100

Paralleling implant gauge




Included in the kit.

Reference NPG 16 100

Gauges and paralleling pin


Gingival and drilling site depth gauge



Not included in the kit.
Drilling depth and gingival height indications.

Reference NJP 22 280


Depth gauge Ø 2.2



Included in the kit.

Reference NJP 22

Naturactis depth gauges



Not included in the kit.

Diameters	References
2.2 - 2.8	NJP 22 28
2.8 - 3.3	NJP 28 33
3.3 - 3.8	NJP 33 38
3.8 - 4.3	NJP 38 43
4.3 - 4.8	NJP 43 48

Paralleling pin Ø 1.8 - 2.2



Included in the kit.

Reference NAP 15 22 18

For common surgical instruments see p 66-69
For common prosthetics instruments see p 174-177


SURGICAL KIT



 Download protocols



Naturactis / Naturall+ common surgical kit



This kit includes all the instruments necessary to place all lengths and all diameters of Naturactis and Naturall+ implants, in all bone densities.


- Simple and compact.
- One kit for 2 implant systems.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use and protocol table.
- Colour coding of plugs according to implant diameter.

Reference NCPT 00

Universal stops kit see p 67

Contents of the kit reference NCPT 00

1	Point drills	\varnothing 1.5 - 2.2	CFP 15 22 50
		\varnothing 2.2	NFP 22 180
2	Initial drills \varnothing 2.2	length 6 mm	AFI 22 060
		length 8 mm	AFI 22 080
		length 10 mm	AFI 22 100
		length 12 mm	AFI 22 120
		length 14 mm	AFI 22 140
		length 16 mm	AFI 22 160
		length 18 mm	AFI 22 180
3	Staged drills \varnothing 2.2 - 2.8	short	NFP 22 28 140
		long	NFP 22 28 180
4	Staged drills \varnothing 2.8 - 3.3	short	NFP 28 33 140
		long	NFP 28 33 180
5	Staged drills \varnothing 3.3 - 3.8	short	NFP 33 38 140
		long	NFP 33 38 180
6	Staged drills \varnothing 3.8 - 4.3	short	NFP 38 43 140
		long	NFP 38 43 180
7	Staged drills \varnothing 4.3 - 4.8	short	NFP 43 48 140
		long	NFP 43 48 180
8	Naturall+ cortical drills	\varnothing 3.5	NFE 35 22
		\varnothing 4	NFE 40 26
		\varnothing 4.5	NFE 45 30
		\varnothing 5	NFE 50 34
9	Naturall+ hard bone drills	\varnothing 3.5	NFC 35 DO
		\varnothing 4	NFC 40 DO
		\varnothing 4.5	NFC 45 DO
		\varnothing 5	NFC 50 DO
10	Naturactis taps	\varnothing 3.5	NTP 35 100
		\varnothing 4	NTP 40 100
		\varnothing 4.5	NTP 45 100
		\varnothing 5	NTP 50 100
Instruments	Depth gauge	\varnothing 2.2	NJP 22
	Paralleling pin	\varnothing 1.8 - 2.2	NAP 15 22 18
	Implants paralleling gauge		NPG 16 100
	Direct implant keys	short	CCP 35 20
		medium	CCP 35 30
		long	CCP 35 40
	Direct implant mandrels	short	CMP 35 20
		medium	CMP 35 30
	External hexagonal keys	short	CCL HE 12 18
		medium	CCL HE 12 22
		long	CCL HE 12 30
	External hexagonal mandrels	short	CMA HE 12 22
		long	CMA HE 12 26
	Extension mandrel		CRM 340
Click wrench		CCC 120	




1

Surgery

Cylindrical implant systems



- 30 Natea+
- 32 Aesthetica+²
- 34 Instruments & common surgical kit
- 42 Uneva+
- 44 Instruments & surgical kit





Natea+

MULTIPURPOSE IMPLANT




-
- Internal hexagonal conical connection
 - the same for all implant diameters
 - compatible with Naturactis+ and Naturall+ implants

-
- Cylindrical implant
 - Bone level

INDICATIONS


- Mandibular arch
- All bone densities, especially high bone densities

Technical characteristics



A	Ø 4.1 - 4.2 - 4.9 - 6.2 mm
B	Smooth neck 0.4 mm
C	Thread 0.3 mm
D	Microthread 3 mm
E	Real screw thread 1.6 mm
F	Thread 0.8 mm
G	Ø 3.6 - 4.1 - 4.8 - 6 mm
H	Sandblasted and etched length

Implant references

 Medical grade titanium.

 Implants are delivered with a cover screw.

Length H	Ø 3.6 mm	Ø 4.1 mm	Ø 4.8 mm	Ø 6 mm
6 mm	-	NIDP 41 42 060	NIDP 48 49 060	NIDP 60 62 060
8 mm	NIDP 36 37 080	NIDP 41 42 080	NIDP 48 49 080	NIDP 60 62 080
10 mm	NIDP 36 37 100	NIDP 41 42 100	NIDP 48 49 100	NIDP 60 62 100
12 mm	NIDP 36 37 120	NIDP 41 42 120	NIDP 48 49 120	NIDP 60 62 120
14 mm	NIDP 36 37 140	NIDP 41 42 140	NIDP 48 49 140	-



Aesthetica+²

TISSUE LEVEL IMPLANT



- Internal octagonal conical connection
- 3 emergence profiles:




- Cylindrical implant
- Tissue level
- One surgical time

INDICATIONS

- Posterior region
- Restoration with wide emergence profile on resorbed ridges
- All bone densities

Technical characteristics



A	Ø 4.2 - 4.8 - 6.5 mm
B	Smooth and etched neck 1.3 mm
C	Thread 0.3 mm
D	Microthread 2.3 mm
E	Real screw thread 1.6 mm
F	Thread 0.8 mm
G	Ø 3.6 - 4.1 - 4.8 mm
H	Sandblasted and etched length

Implant references

Medical grade titanium.


Implants are delivered with a 3 mm in height healing abutment.

Length H	Ø 3.6 mm		Ø 4.1 mm		Ø 4.8 mm	
	neck Ø 4.2 mm	neck Ø 4.8 mm	neck Ø 4.2 mm	neck Ø 4.8 mm	neck Ø 4.8 mm	neck Ø 6.5 mm
6 mm	-	-	AIEP 41 42 060	AIEP 41 48 060	AIEP 48 48 060	AIEP 48 65 060
8 mm	AIEP 36 42 080	AIEP 36 48 080	AIEP 41 42 080	AIEP 41 48 080	AIEP 48 48 080	AIEP 48 65 080
10 mm	AIEP 36 42 100	AIEP 36 48 100	AIEP 41 42 100	AIEP 41 48 100	AIEP 48 48 100	AIEP 48 65 100
12 mm	AIEP 36 42 120	AIEP 36 48 120	AIEP 41 42 120	AIEP 41 48 120	AIEP 48 48 120	AIEP 48 65 120
14 mm	AIEP 36 42 140	AIEP 36 48 140	AIEP 41 42 140	AIEP 41 48 140	AIEP 48 48 140	AIEP 48 65 140

INSTRUMENTS

Drills


Point drill Ø 1.5 - 2.2



The upper part of the point drill has the same diameter 2.2 mm of the following drill to prepare its insertion.
The 90° cutting edge makes a clear and precise mark on the bone crest, even when the crest is thin.
Included in the kit.

Reference CFP 15 22 50


Initial drills Ø 2.2



Included in the kit.

Lengths	References
6	AFI 22 060
8	AFI 22 080
10	AFI 22 100
12	AFI 22 120
14	AFI 22 140


Staged drills



Included in the kit.

Diameters	Lengths	References
Ø 2.2 - 2.8	14	NFP 22 28 140
	18	NFP 22 28 180
Ø 2.8 - 3.3	14	NFP 28 33 140
	18	NFP 28 33 180
Ø 3.3 - 3.8	14	NFP 33 38 140
	18	NFP 33 38 180
Ø 3.8 - 4.3	14	NFP 38 43 140
	18	NFP 38 43 180

Straight drills for Natea+ Ø 6



Included in the supplementary surgical kit Natea+ Ø 6.

Diameters	References
4.8	AFD 48 120
5.2	AFD 52 120
5.4	AFJ 54 120
5.7	AFD 57 120

Drills - continuation


Cortical drills



Included in the kit and in the supplementary surgical kit Natea+ Ø 6.

Implants	References
Natea+ Ø 3.6 / 3.7	ANFP 36 30
Aesthetica+ ² Ø 3.6 / 4.2	ANFP 42 30
Natea+ / Aesthetica+ ² Ø 4.2	ANFP 42 35
Natea+ / Aesthetica+ ² Ø 4.8 / 4.8	ANFP 48 42
Aesthetica+ ² Ø 4.8 / 6.5	ANFP 48 43
Natea+ Ø 6	ANFP 60 52


Taps



Included in the kit and in the supplementary surgical kit Natea+ Ø 6.

Diameters	References
3.6	ATB 36 126
4.1	ATR 41 126
4.8	ATV 48 126
6	ATJ 60 120

Implants paralleling pin




Included in the kit.

Implants	References
Natea+	NPG 16 100
Aesthetica+ ²	APG 20 100

Gauges and paralleling pins - continuation


Gingival and drilling site depth gauges



Not included in the kit.
Drilling depth and gingival height indications.

Implants	References
Natea+	NJP 22 280
Aesthetica+ ²	AJP 22 280

Depth gauge Ø 2.2



Included in the kit.

Reference NJP 22


Paralleling pin Ø 1.8 - 2.2



Included in the kit.

Reference NAP 15 22 18

Paralleling pins



Not included in the kit.

Shoulder Ø a	References
4.8	AAP 48 22 28
6.5	AAP 65 22 28

For common surgical instruments see p 66-69

For common prosthetics instruments see p 174-177


SURGICAL KITS



 Download protocols



Natea+ / Aesthetica+2 Ø 3.6 - 4.1 - 4.8 common surgical kit



This kit includes all the instruments necessary to place all lengths of Natea+ and Aesthetica+2 Ø 3.6 - 4.1 - 4.8. implants, in all bone densities.

- Simple and compact.
- One kit for 2 implant systems.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use and protocol table.
- Colour coding of plugs according to implant diameter.

Reference ANCP T 00

Universal stops kit see p 67

Contents of the kit reference ANCPT 00


1	Point drill Ø 1.5 - 2.2		CFP 15 22 50
2	Initial drills Ø 2.2	<i>length 6 mm</i>	AFI 22 060
		<i>length 8 mm</i>	AFI 22 080
		<i>length 10 mm</i>	AFI 22 100
		<i>length 12 mm</i>	AFI 22 120
		<i>length 14 mm</i>	AFI 22 140
3	Staged drills 2.2 - 2.8	<i>short</i>	NFP 22 28 140
		<i>long</i>	NFP 22 28 180
4	Staged drills Ø 2.8 - 3.3	<i>short</i>	NFP 28 33 140
		<i>long</i>	NFP 28 33 180
5	Staged drills Ø 3.3 - 3.8	<i>short</i>	NFP 33 38 140
		<i>long</i>	NFP 33 38 180
6	Staged drills Ø 3.8 - 4.3	<i>short</i>	NFP 38 43 140
		<i>long</i>	NFP 38 43 180
7	Cortical drills	<i>Natea+ - Ø 3.6 / 3.7</i>	ANFP 36 30
		<i>Aesthetica+² - Ø 3.6 / 4.2</i>	ANFP 42 30
		<i>Natea+ / Aesthetica+² - Ø 4.1</i>	ANFP 42 35
		<i>Natea+ / Aesthetica+² - Ø 4.8 / 4.8</i>	ANFP 48 42
		<i>Aesthetica+² Ø 4.8 / 6.5</i>	ANFP 48 43
8	Taps	<i>Ø 3.6</i>	ATB 36 126
		<i>Ø 4.1</i>	ATR 41 126
		<i>Ø 4.8</i>	ATV 48 126
Instruments	Depth gauge	<i>Ø 2.2</i>	NJP 22
	Paralleling pin	<i>Ø 1.8 - 2.2</i>	NAP 15 22 18
	Implants paralleling gauges	<i>Aesthetica+²</i>	APG 20 100
		<i>Natea+</i>	NPG 16 100
	Direct implant keys	<i>short - Natea+</i>	CCP 35 20
		<i>medium - Natea+</i>	CCP 35 30
		<i>long - Natea+</i>	CCP 35 40
	Direct implant mandrels	<i>short - Natea+</i>	CMP 35 20
		<i>medium - Natea+</i>	CMP 35 30
	Direct implant keys	<i>short - Aesthetica+²</i>	CCP 42 20
		<i>medium - Aesthetica+²</i>	CCP 42 30
		<i>long - Aesthetica+²</i>	CCP 42 40
	Direct implant mandrels	<i>short - Aesthetica+²</i>	CMP 42 20
		<i>medium - Aesthetica+²</i>	CMP 42 30
	External hexagonal keys	<i>short</i>	CCL HE 12 18
		<i>medium</i>	CCL HE 12 22
		<i>long</i>	CCL HE 12 30
	External hexagonal mandrels	<i>short</i>	CMA HE 12 22
<i>long</i>		CMA HE 12 26	
Extension mandrel		CRM 340	
Click wrench		CCC 120	



 Download protocol



Natea+ Ø 6 supplementary surgical kit



This kit includes the instruments complementary to the kit ANCPT 00 for surgery and the management of all bone densities and all lengths of Natea+ Ø 6 mm implants.

- Simple and compact.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use and protocol table.

Reference ANCPT 60

Universal stops kit see p 67

Contents of the supplementary surgical kit reference ANCPT 60

1	<i>Cylindrical drills</i>	\varnothing 4.8	AFD 48 120
		\varnothing 5.2	AFD 52 120
		\varnothing 5.4	AFJ 54 120
2	<i>Terminal drill</i>		AFD 57 120
3	<i>Tap</i>		ATJ 60 120
4	<i>Cortical drill</i>		ANFP 60 52



Uneva+

EXTERNAL CONNECTION IMPLANT




-
- External hexagonal connection

-
- Cylindrical implant
 - Bone level

INDICATIONS

- Multipurpose implant
- Suitable for high bone densities
- Direct bars on implants

Technical characteristics



A	Ø 4.1 - 4.2 - 4.9 mm
B	Smooth neck 0.4 to 0.8 mm according to the implant diameter
C	Thread 0.3 mm
D	Microthread 2.9 mm
E	Real screw thread 1.6 mm
F	Thread 0.8 mm
G	Ø 3.6 - 4.1 - 4.8 mm
H	Sandblasted and etched length

Implant references

 Medical grade titanium.


 Implants are delivered with a cover screw.

Length H	Ø 3.6 mm	Ø 4.1 mm	Ø 4.8 mm
6 mm	-	-	UHDP 48 49 060
8 mm	UHDP 36 41 080	UHDP 41 41 080	UHDP 48 49 080
10 mm	UHDP 36 41 100	UHDP 41 41 100	UHDP 48 49 100
12 mm	UHDP 36 41 120	UHDP 41 41 120	UHDP 48 49 120
14 mm	UHDP 36 41 140	UHDP 41 41 140	UHDP 48 49 140

INSTRUMENTS

Drills


Point drill Ø 1.5 - 2.2



The upper part of the point drill has the same diameter 2.2 mm of the following drill to prepare its insertion.
The 90° cutting edge makes a clear and precise mark on the bone crest, even when the crest is thin.
Included in the kit.

Reference CFP 15 22 50


Initial drills Ø 2.2



Included in the kit.

Lengths	References
6	AFI 22 060
8	AFI 22 080
10	AFI 22 100
12	AFI 22 120
14	AFI 22 140


Staged drills



Included in the kit.

Diameters	Lengths	References
Ø 2.2 - 2.8	14	NFP 22 28 140
	18	NFP 22 28 180
Ø 2.8 - 3.3	14	NFP 28 33 140
	18	NFP 28 33 180
Ø 3.3 - 3.8	14	NFP 33 38 140
	18	NFP 33 38 180
Ø 3.8 - 4.3	14	NFP 38 43 140
	18	NFP 38 43 180


Cortical drills



Included in the kit.

Diameters	References
3.6	ANFP 41 30
4.1	ANFP 42 35
4.8	ANFP 48 42


Taps



Included in the kit.

Diameters	References
3.6	ATB 36 126
4.1	ATR 41 126
4.8	ATV 48 126

Implants paralleling pin




Included in the kit.

Reference UPG 20 100

Gauges and paralleling pins

Gingival and drilling site depth gauge




Not included in the kit.
Drilling depth and gingival height indications.

Reference UJP 22 280

Gauges and paralleling pins - continuation

Depth gauge Ø 2.2



Included in the kit.

Reference NJP 22

Paralleling pin Ø 1.8 - 2.2



Included in the kit.

Reference NAP 15 22 18

Paralleling pin Ø 4.8



Not included in the kit.

Reference AAP 48 22 28

For common surgical instruments see p 66-69
For common prosthetics instruments see p 174-177


SURGICAL KIT



 Download protocol



Uneva+ surgical kit



This kit includes all the instruments necessary to place all lengths and all diameters of Uneva+ implants, in all bone densities.



- Simple and compact.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use and protocol table.
- Colour coding of plugs according to implant diameter.

Reference UCPT 00

Universal stops kit see p 67

Contents of the kit reference UCPT 00


1	Point drill Ø 1.5 - 2.2		CFP 15 22 50
2	Initial drills Ø 2.2	<i>length 6 mm</i>	AFI 22 060
		<i>length 8 mm</i>	AFI 22 080
		<i>length 10 mm</i>	AFI 22 100
		<i>length 12 mm</i>	AFI 22 120
		<i>length 14 mm</i>	AFI 22 140
3	Staged drills 2.2 - 2.8	<i>short</i>	NFP 22 28 140
		<i>long</i>	NFP 22 28 180
4	Staged drills Ø 2.8 - 3.3	<i>short</i>	NFP 28 33 140
		<i>long</i>	NFP 28 33 180
5	Staged drills Ø 3.3 - 3.8	<i>short</i>	NFP 33 38 140
		<i>long</i>	NFP 33 38 180
6	Staged drills Ø 3.8 - 4.3	<i>short</i>	NFP 38 43 140
		<i>long</i>	NFP 38 43 180
7	Cortical drills	Ø 3.6	ANFP 41 30
		Ø 4.1	ANFP 42 35
		Ø 4.8	ANFP 48 42
8	Taps	Ø 3.6	ATB 36 126
		Ø 4.1	ATR 41 126
		Ø 4.8	ATV 48 126
Instruments	Depth gauge	Ø 2.2	NJP 22
	Paralleling pin	Ø 1.8 - 2.2	NAP 15 22 18
	Implants paralleling gauge		UPG 20 100
	Direct implant keys	<i>short</i>	CCP 24 20
		<i>medium</i>	CCP 24 30
		<i>long</i>	CCP 24 40
	Direct implant mandrels	<i>short</i>	CMP 24 20
		<i>medium</i>	CMP 24 30
	External hexagonal keys	<i>short</i>	CCL HE 12 18
		<i>medium</i>	CCL HE 12 22
		<i>long</i>	CCL HE 12 30
	External hexagonal mandrels	<i>short</i>	CMA HE 12 22
		<i>long</i>	CMA HE 12 26
	Extension mandrel		CRM 340
Click wrench		CCC 120	



1

Surgery

Narrow implant systems



- 52 Naturactis Ø 3
- 54 Naturall+ Ø 3
- 56 Instruments & common surgical kit
- 60 Obi Ø 2.7
- 62 Instruments & surgical kit




Naturactis Ø3 NARROW IMPLANT

1

Naturactis Ø3

SURGERY - Implant systems - Narrow implants




-
- Internal hexagonal conical connection compatible with the Naturall+ Ø 3 implant

-
- Tapered implant
 - Bone level

INDICATIONS

- Restorations of inferior, central mandibular and lateral incisors
- Reduced mesiodistal spaces
- Thin ridges in the anterior region
- Post-extraction surgery
- Immediate loading

Technical characteristics



A	Ø 3
B	Supracrestal smooth neck 0.4 mm
C	Thread 0.4 mm
D	Microthread 3 mm
E	Real screw thread 2.4 mm
F	Thread 1.2 mm
G	Ø 1.5 mm
H	Sandblasted and etched length

Implant references

 TA6V ELI medical grade.

 Implants are delivered with a cover screw.

Length H	Ø 3 mm
8 mm	NIP 30 080
10 mm	NIP 30 100
12 mm	NIP 30 120
14 mm	NIP 30 140




Naturall+ Ø3

NARROW IMPLANT

1

Naturall+ Ø3

SURGERY - Implant systems - Narrow implants




-
- Internal hexagonal conical connection compatible with the Naturactis Ø 3 implant

-
- Tapered implant
 - Bone level

INDICATIONS

- Restorations of inferior, central mandibular or lateral incisors
- Reduced mesiodistal spaces
- Thin ridges in the anterior region
- Post-extraction surgery
- Immediate loading

Technical characteristics



A	Ø 3 mm
B	Supracrestal smooth neck 0.4 mm
C	Thread 0.3 mm
D	Microthread 3 mm
E	Real screw thread 1.8 mm
F	Thread 0.9 mm
G	Ø 2.2
H	Sandblasted and etched length

Implant references

 TA6V ELI medical grade.


 Implants are delivered with a cover screw.

Length H	Ø 3 mm
8 mm	NICP 30 080
10 mm	NICP 30 100
12 mm	NICP 30 120
14 mm	NICP 30 140

INSTRUMENTS

Drills


Point drill Ø 1.5 - 2.2



The upper part of the point drill has the same diameter 2.2 mm of the following drill to prepare its insertion.
The 90° cutting edge makes a clear and precise mark on the bone crest, even when the crest is thin.
Included in the kit.

Reference CFP 15 22 50


Initial drill Ø 1.8



Included in the kit.

Reference NFI 18 150


Drill Ø 2.2



Included in the kit.

Reference AFI 22 140


Short staged drill Ø 2.2 - 2.8



Included in the kit.

Reference NFP 22 28 140


Cortical drill



Included in the kit.

Reference NFE 30 18

Naturall+ hard bone drill Ø 3




Included in the kit.

Reference NFC 30 DO

Depth gauge and paralleling pin

Depth gauge Ø 1.8



Not included in the kit.

Reference NJP 18 250

Paralleling pin Ø 1.8 - 2.2




Included in the kit.

Reference NAP 15 22 18

For common surgical instruments see p 66-69
 For common prosthetics instruments see p 174-177


SURGICAL KIT



 Download protocols



Naturactis Ø 3 / Naturall+ Ø 3 common surgical kit



This kit includes all the instruments necessary to place all lengths of Naturactis Ø 3 and Naturall+ Ø 3 implants, in all bone densities.

- Simple and compact.
- One kit for 2 implant systems.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use and protocol table.

Reference NCPT 30

Universal stops kit see p 67


Contents of the kit reference NCPT 30

1	<i>Point drill</i>	Ø 1.5 - 2.2	CFP 15 22 50
2	<i>Initial drill</i>	Ø 1.8	NFI 18 150
3	<i>Cortical drill</i>	Ø 3	NFE 30 18
4	<i>Drill</i>	Ø 2.2	AFI 22 140
5	<i>Naturall + hard bone drill</i>	Ø 3	NFC 30 DO
	<i>Short staged drill</i>	Ø 2.2 - 2.8	NFP 22 28 140
<i>Instruments</i>	<i>Depth gauge</i>	Ø 1.8	NJP 18 250
	<i>Paralleling pin</i>	Ø 1.8 - 2.2	NAP 15 22 18
	<i>Direct implant key</i>	long	CCP 30 40
	<i>Direct implant mandrel</i>	long	CMP 30 30
	<i>External hexagonal key</i>	long	CCL HE 12 30
	<i>External hexagonal mandrel</i>	long	CMA HE 12 26
	<i>Extension mandrel</i>		CRM 340
	<i>Click wrench</i>		CCC 120



Obi Ø2.7

MINI-IMPLANT




- Simple surgical protocol = only one drill required
- Universal ball for O-Ring Ø 2.25 mm

INDICATIONS

- Stabilization of removable prosthesis


Technical characteristics



A	Ø 3.5 mm
B	4 mm
C	0.5 mm
D	2.8 mm
E	Ø 2.7 mm
F	Real screw thread 1.6 mm
G	Thread 0.8 mm
H	Ø 2.35 mm
I	Sandblasted and etched length

Implant references

 TA6V ELI medical grade.


 Implants are delivered with a 60 shores O-Ring attachment.

Length I	Ø 2.7
9 mm	OIC 27 68 090
11 mm	OIC 27 68 110
13 mm	OIC 27 68 130
15 mm	OIC 27 68 150

INSTRUMENTS

Gingival cutter and drills


Gingival cutter



Included in the kit.

Reference ODG 20 35

Point drill Ø 1.5 - 2.2




The upper part of the point drill prepare the insertion of the following drill.

The 90° cutting edge makes a clear and precise mark on the bone crest, even when the crest is thin.

Included in the kit.

Reference CFP 15 22 50

Drill Ø 2




Included in the kit.

Reference OFI 20 150

Depth gauge and paralleling pin

Depth gauge Ø 1.8



Included in the kit.

Reference NJP 18 250

Paralleling pin Ø 1.8 - 2.2



Included in the kit.

Reference NAP 15 22 18

For common surgical instruments see p 66-69
For common prosthetics instruments see p 174-177

SURGICAL KIT

1

Instruments and surgical kit


SURGERY - Implant systems - Narrow implants - Obi Ø 2.7



i Download protocol



Obi Ø 2.7 surgical kit



This kit includes all the instruments necessary to place all lengths of Obi Ø 2.7 implants, in all bone densities.

- Simple and compact.
- Reduced format for increased space on a sterile field and in an autoclave.
- Legibility of sequences thanks to the layout of the instruments in order of use.

Reference OICK 27 XX 00

Contents of the kit reference OICK 27 XX 00

1	Gingival cutter		ODG 20 35
2	Point drill	Ø 1.5 - 2.2	CFP 15 22 50
3	Drill	Ø 2	OFI 20 150
Instruments	Depth gauge	Ø 1.8	NJP 18 250
	Paralleling pin		NAP 15 22 18
	Internal hexagonal mandrel		CMO HI 25 26
	Internal hexagonal keys	short	CCL HI 25 18
		long	CCL HI 25 26
	Click wrench		CCC 120



1

Surgery

Peripheral products




- 66 Common instruments
- 70 Extraction kits
- 72 Bone filling Macrobone®



COMMON INSTRUMENTS

Surgical torque wrench



For placing implants in the surgical site.
 With progressive torque adjustment every 5 Ncm, within the recommended range of 35 Ncm to 75 Ncm. The torque wrench can be used instead of a surgical motor when the surgeon prefers to insert the implant by hand or for motors that do not display the torque level.
 Allows accurate measurement of primary implant stability.

- Made of surgical stainless steel.
- Detachable for easier cleaning.
- Different torque adjustments available: 35, 40, 45, 50, 55, 60, 65, 70 & 75 Ncm.
- Automatic stop: tightening stops automatically when the pre-set value is reached.




Download the cleaning instructions



Reference CCC 70

Surgical sequencer



This allows you to lay out the implants and instruments in the right order for a specific surgical procedure.
 DELIVERED EMPTY.

Example of use

Reference CSC 7 20

Stops

Stops kit





Universal kit: can be used with all implant ranges (only for the 2 level drills).

- Gripping of stops directly with contra-angle.
- Colour coding for easy identification of stops according to the implant to be fitted.
- 28 stops for short and long drills included in the kit.
- Kit can be sterilised in an autoclave.

Reference KBU 00
Empty delivered reference KBU 00V

Stops



Sold individually.
 Colour coding.

For drills diameters	Lengths mm	References
2.2 - 3	2	CBU 33 02
	4	CBU 33 04
	6	CBU 33 06
	8	CBU 33 08
	10	CBU 33 10
	12	CBU 33 12
3.3 - 3.5	14	CBU 33 14
	2	CBU 40 02
	4	CBU 40 04
	6	CBU 40 06
	8	CBU 40 08
	10	CBU 40 10
3.8 - 4.5	12	CBU 40 12
	14	CBU 40 14
	2	CBU 47 02
	4	CBU 47 04
	6	CBU 47 06
	8	CBU 47 08
4.8 - 5.7	10	CBU 47 10
	12	CBU 47 12
	14	CBU 47 14
	2	CBU 60 02
	4	CBU 60 04
	6	CBU 60 06
	8	CBU 60 08
	10	CBU 60 10
	12	CBU 60 12
	14	CBU 60 14

Stop removal tool




Included in the kit.

Reference ACB 36 48


Keys and mandrels

Direct implant keys



	1. Naturactis Ø 3 Naturall+ Ø 3	2. Naturactis Naturall+ Natea+	3. Aesthetica+ ²	4. Uneva+
Short - L 10	CCP 30 20	CCP 35 20	CCP 42 20	CCP 24 20
Medium - L 15	CCP 30 30	CCP 35 30	CCP 42 30	CCP 24 30
Long - L 20	CCP 30 40	CCP 35 40	CCP 42 40	CCP 24 40

Direct implant mandrels



	1. Naturactis Ø 3 Naturall+ Ø 3	2. Naturactis Naturall+ Natea+	3. Aesthetica+ ²	4. Uneva+
Short - L 30	CMP 30 20	CMP 35 20	CMP 42 20	CMP 24 20
Long - L 35	CMP 30 30	CMP 35 30	CMP 42 30	CMP 24 30

Square tap key



For tap manual use.
Allows tap to be used with click wrench and torque wrench.

Reference CEC 40

Mandrel extension



Use as extension for all drills and instruments that can be used with a contra-angle.

Reference CRM 340

Click wrench




Detachable.

Reference CCC 120


Keys and mandrels - continuation

Screwing handle for Naturactis - Naturall+ - Natea+ implants



Reference
CAM 70 100


Drilling guides



Point of reference for the drill axis.
Indication of the minimum distance between 2 implants with a vertical line every 7 mm, reproducing average interdental space.
To be used after initial drilling.

	Ø	References
For Obi	2	CGF 20 650
Universal	2.2	CGF 22 650


Tablets for testing titanium allergy



To carry out sensitivity tests on patients by an allergist.

	References
Grade IV titanium	ETK TIG 4
Grade V titanium	ETK TIG 5

Sinus parachute screw



This screw is placed on implants inserted into the maxilla close to the sinus where stability may be sufficient to prevent them in the long term from sliding into sinus cavities and potentially perforating the sinus membrane.

	Reference
Naturall+ Ø 3.5 - 4 - 4.5	NVS 35 50
Natea+ Ø 3.6 - 4.1	

EXTRACTION KITS

A complete solution for mechanical complications encountered in implantology.



1 - Removal of blocked abutments in an implant

Abutments can be removed without interfering with the implant thanks to an abutment extractor.

2 - Removal of a broken screw in an implant

- Thanks to a drilling guide and a carbide drill with a high cutting ability of titanium, a small opening is made in the remaining part of the screw.


- A reverse drill is then inserted into the opening in the broken screw to allow unscrewing to take place.

- A tap can be used to retouch the threading of the implant if it is damaged.


3 - Removal of an implant

- The implant can be unscrewed using an implant extractor.

- The bone around the implant can be trepanned using trephines suitable for the diameters of the different implants.



1



2









3

Kits and options/accessories references

Ranges	Naturactis / Naturall+ / Natea+	Aesthetica+ ²	Uneva+
Kits references	KDR 3N	KDR AEST	KDR U
Options/accessories references	Aesthetica+ ² : KIR AEST	Naturactis / Naturall+ / Natea+: KIR 3N	Naturactis / Naturall+ / Natea+: KIR 3N
	Uneva+: KIR U	Uneva+: KIR U	Aesthetica+ ² : KIR AEST

Photo non contractuelle

Kits and options/accessories contents

	References	Kits			Options/ accessories		
		KDR 3N	KDR AEST	KDR U	KIR 3N	KIR AEST	KIR U
	M1.6 abutment extractor CEP 16 175	X			X		
	M2 abutment extractor CEP 20 175		X	X		X	X
	Reverse drill Ø 1.2 CFK TG 12	X	X	X			
	Naturactis / Naturall+ / Natea+ guide CGK 3N	X			X		
	Aesthetica+² guide CGK AEST		X			X	
	Uneva+ guide CGK U			X			X
	M1.6 tap CTK TA 16	X			X		
	M2 tap CTK TA 20		X	X		X	X
	Screw extractor CCK QV 15 150	X	X	X			
	Implant extractor CCK QI 15 150	X	X	X			
	Trepine Ø 3.8 CMK TR 38	X	X	X			
	Trepine Ø 4.6 CMK TR 46	X	X	X			
	Trepine Ø 5.3 CMK TR 53	X	X	X			
	Trepine Ø 6.7 CMK TR 67	X	X	X			

Extraction kits

SURGERY - Peripheral products

BONE FILLING


Macrobone

Ultra-porous synthetic bone (close to 90% porosity) designed for filling bone defects in oral and maxillofacial surgery.

Composition: pure-phase beta-tricalcium phosphate


β -TCP = a bioactive calcium phosphate salt that is more soluble in a biological environment than hydroxyapatite and very close to human bone mineral.

- **Biocompatible** = osseointegrated then resorbed.
- **The material is totally resorbed** and replaced by newly-formed bone in 6-7 months.
- **Completely synthetic** = no immune or infection risks.
- **Radio-opaque** = radiological monitoring of the cavity filling and graft integration.



Porosity close to 90%


- **Macroporosity** = large pores (0.2–0.5 mm) that are interconnected with each other, allowing the bone to penetrate into the heart of the material (osseointegration by osteoinduction).
- **Microporosity.**
- **Presence of biomaterials is limited to 10%** compared with its total volume due to high porosity.
 - Faster transformation.
 - Easier cell progression.
- **Asymmetrical pellet shape** = no problems filling irregularly shaped cavities.



High porosity of 90%. French Institute of Health and Medical Research (INSERM) - Angers (Professor D. Chappard)

Histology and clinical experience

- Study carried out by the French Institute of Health and Medical Research (Institut national de la santé et de la recherche médicale - INSERM)
Method => Bone core samples taken 3 months after sinus lift with Macrobone®.
Result => All showed the same bone remodeling appearance. The image opposite shows the perfect placement of the bone on the surface of the biomaterial, signifying its osseointegration.
- Macrobone® benefits from **8 years' clinical experience**.



Coloration: Goldner's trichrome
- - - - - : limit of the implanted pellet

Direct positioning of the bone on the surface of the bone substitute.


Range and indications



Pellet size:
150-500 µm (0,15 to 0,5 mm)
Indications: periodontology.



Pellet size:
500-1000 µm (0,5 to 1 mm)
Indications: socket filling; filling of medium-sized bone defects.




Pellet size:
1000-2000 µm (1 to 2 mm)
Indications: sinus lift.

Packaging and references




i Packaged in a box in 2 individual sterile wells (double sterile packaging).

	Contents	
Pellet size	2 wells x 0,5cc	2 wells x 1cc
150-500 µm	BON 15 050 052	BON 15 050 102
500-1000 µm	BON 50 100 052	BON 50 100 102
1000-2000 µm	BON 100 200 052	BON 100 200 102



Part 2



Healing




- 76 Titanium healing abutments and cover screws
- 80 Anatomical healing (Profile Designer iphysio®)

TITANIUM HEALING ABUTMENTS & COVER SCREWS

Choice of healing abutment

The healing abutment shapes the future emergence profile of the implant when the gingiva is healing.

Choosing the correct healing abutment depends on:
 A - the desired emergence of the permanent abutment,
 B - the depth of the prosthetic/abutment junction.




Healing abutments have a diameter that is slightly greater (0.4 mm) than the permanent abutment:

- to ensure easier and less painful insertion of impression copings and abutments (avoiding the need to anaesthesia),
- to avoid entrapping the gingiva and improving patient comfort,
- to speed up the procedure.

NATURACTIS - NATURALL+ - NATEA+

Healing abutments



Healing abutment and final abutment height correspondence



Titanium.



* Do not use with Naturactis implant.

Prosthetic profiles	Healing abutments references	Supra-implant heights D
	NCI 36 23	0.5
	NCI 36 34	1.5
	NCI 36 45	2.5
	NCI 36 56	3.5
	NCI 36 67	4.5
	NCI 46 23	0.5
	NCI 46 34	1.5
	NCI 46 45	2.5
	NCI 46 56	3.5
	NCI 46 67	4.5
	NCI 52 23	0.5
	NCI 52 34	1.5
	NCI 52 45	2.5
	NCI 52 56	3.5
	NCI 52 67	4.5
	NCI 60 34*	1.5
	NCI 60 45	2.5
	NCI 60 56	3.5
	NCI 60 67	4.5

Cover screw



Titanium.




Delivered with the implant.

Reference NVC 35 01

NATURACTIS Ø 3 - NATURALL+ Ø 3

Healing abutments



Healing abutment and final abutment height correspondence



Titanium.



Prosthetic profile	Healing abutments references	Supra-implant height D
	NCI 30 23	0.5
	NCI 30 45	2.5
	NCI 30 67	4.5

NATURACTIS Ø 3 - NATURALL+ Ø 3 - CONTINUATION

Cover screw



Titanium.




Delivered with the implant.



Reference NVC 31

AESTHETICA+²

Healing abutments



Healing abutments (height = 3 mm) are supplied with the implant.



Titanium.



Prosthetic profiles	Supra-implant heights	References
	1.5	ACI 42 47 15
	3	ACI 42 47 30
	4.5	ACI 42 47 45
	1.5	ACI 48 55 15
	3	ACI 48 55 30
	4.5	ACI 48 55 45
	2	ACI 65 72 20
	3	ACI 65 72 30
	4.5	ACI 65 72 45

Cover screws




Titanium.



Implant diameters	References
4.2	ACI 48 35 00
4.8	
6.5	ACI 65 43 00

UNEVA+




Healing abutments



Healing abutment and final abutment height correspondence

 Titanium.

 10Ncm

Implant diameters	Prosthetic profiles	Supra-implant heights	References
3.6 - 4.1		1	UCI 37 46 25
		2	UCI 37 46 35
		3	UCI 37 46 45
		1	UCI 37 51 25
		2	UCI 37 51 35
		3	UCI 37 51 45
4.8		1	UCI 47 58 25
		2	UCI 47 58 35
		3	UCI 47 58 45

Cover screws




 Titanium.

 Delivered with the implant.

 10Ncm


Implant diameters	References
3.6 - 4.1	UCO 37
4.8	UCO 47



PROFILE DESIGNER



iphysio®
Profile Designer
Healing
+
Impression





Esthetibase®
interface
Support
iphysio®
+
Prosthesis support



2 in 1

One piece for the healing and impression
= no destruction of the mucosa attachment,
= time saving.



1. Soft tissue healing according to the anatomy of the tooth to replace.


2. Classic or digital impression directly on the cap.

Anatomical


9 references to better match the anatomy of the tooth to replace.

3 profiles

- A** Incisors and canines
- B** Premolars
- C** Molars



3 heights




This system is 3Shape® compatible.

Universal

It can be used with most of implant brands thanks to a wide range of compatible titanium Esthetibase® interfaces.

 Download the instructions of use




 For all Esthetibase interfaces available, see p 162-163 and 180-181




References

iphysio® small platform on size S Esthetibase interface



	Shape A	Shape B	Shape C
Height 1 mm	ETK_EB.S.A1	ETK_EB.S.B1	ETK_EB.S.C1
Height 2 mm	ETK_EB.S.A2	ETK_EB.S.B2	ETK_EB.S.C2
Height 4 mm	ETK_EB.S.A4	ETK_EB.S.B4	ETK_EB.S.C4

iphysio® large platform on size L Esthetibase interface



	Shape A	Shape B	Shape C
Height 1 mm	ETK_EB.L.A1	ETK_EB.L.B1	ETK_EB.L.C1
Height 2 mm	ETK_EB.L.A2	ETK_EB.L.B2	ETK_EB.L.C2
Height 4 mm	ETK_EB.L.A4	ETK_EB.L.B4	ETK_EB.L.C4

iphysio® kits




These kits allow you to classify your Profile Designers according to their shapes and their heights.

1 kit for each platform (S and L) which contains:

- 1 scanbody
- 3 test pieces (shapes A - B - C)
- 9 iphysio® Profile Designers

Size	Kits delivered full references	Kits delivered empty references
S	KIP 00 S	KIP 00 SV
L	KIP 00 L	KIP 00 LV



Part 3

etk prosthesis



82 Standard prosthetic components

- 84 Naturactis / Naturall+ / Natea+
- 110 Aesthetica+²
- 126 Uneva+
- 144 Naturactis / Naturall+ Ø 3
- 156 Obi Ø 2.7

160 Customised prosthesis CAD-CAM

- 162 Titanium Esthetibase interfaces
- 164 Scanbody
- 166 CAD-CAM works

168 All in bar[®] system

172 Common instruments

- 174 Prosthetic torque wrench
- 175 Prosthetic kit
- 176 Keys & mandrels



3

etk prosthesis

Standard prosthetic components

.....

Naturactis
Naturall+
Natea+



COMMON PROSTHETIC RANGE FOR NATURACTIS / NATURALL+ / NATEA+

Internal hexagonal conical connection

- Better sealing of the restoration.
- Stability of the implant/component assembly.
- Precision orientation of prosthetic elements.

Single connection for all implant diameters

- The choice of prosthetic platform does not depend on the choice of implant diameter.

Unique connection for 3 implant systems

- Streamlining of stocks of prosthetic components for simplified component management.

A single prosthetic range and 4 prosthetic platforms

- To simplify your treatment plan and facilitate your surgeries and prosthetic constructions.

Connection proven for 10 years

- Proven mechanical resistance.
- Fatigue tests for 5 million cycles according to the ISO 14801 standard.



Height details of prosthetic components have been changed.


The heights indicated are the heights above the implants and are identical for the Naturactis, Naturall+ and Natea+ implants.

4 PROSTHETIC PLATFORMS




1 COMMON
PROSTHETIC RANGE


1 SINGLE
CONNECTION



3 IMPLANT SYSTEMS

PROSTHETIC RECOMMENDATIONS





(1) Within the limit of prosthetic axes divergences and of a reduced number of elements (8 maximum).

CEMENTED RESTORATION ON STRAIGHT AND ANGULATED ABUTMENTS

For single and multi-unit prosthesis




- 4 prosthesis platforms available:




- 5 supra-implant heights: 0.5 / 1.5 / 2.5 / 3.5 / 4.5 mm.
- 3 angles: 7°, 15° and 20°.
- Nitrided abutments for a better aesthetic result.
- Abutments delivered with a titanium fixing screw: the screw is secured thanks to an interior threading of the prosthetic part to avoid sudden collapse.
- Fixing screw with “anti-unscrewing” treatment.
- Laser marking to identify the prosthesis platform and supra-implant height of the abutments.

Letter	Emergence diameter
E	3.6
N	4.6
R	5.2
W	6



Temporization

Direct temporary abutment



Supplied with a translucent white rotational cap (ref. APS CP 36 40) and an opaque white non-rotational cap (ref. APS CO 36 40).
ht = supra-implant height.



Titanium + medical polymer.



Reference NPS PPT 36 40

Temporary trans-screwed rotational abutment



Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment.



For multi-unit restorations.
ht = supra-implant height.



Titanium.

Reference NPS PPTR 38 2

Temporary trans-screwed non-rotational abutment



Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment.



For single restorations.
ht = supra-implant height.




Titanium.

Reference NPS PPT 38 2

Impression

Pop-in impression coping



Supplied with a titanium impression coping screw (ref. NPS VTB 16 156).



Closed tray impression taking.



Titanium.

Reference NPI 37

Pop-up impression coping



Supplied with a titanium impression coping screw (ref. NPS VTB 16 174).



Closed tray impression taking.
Supplied with an impression cap for pick-up technique (ref. APS TCP 36 40).



Titanium + medical polymer.

Reference NPU 35

Impression - continuation

Short pick-up impression coping



Supplied with a titanium laboratory screw (ref. NPS VG 16 200).



Open tray impression taking.



Titanium.

Reference NPE T35

Long pick-up impression coping



Supplied with a titanium laboratory screw (ref. NPS VG 16 250).



Open tray impression taking.



Titanium.

Reference NPE T35 L

Implant analog




Can be cut for use with pins.



Titanium.

Reference NLA H35


Bite registration abutment



Medical polymer.

References	
Sold single	ETK NA 35 PO
Kit of 10	ETK NA 35 PO 10

Scanbody direct on implant



Supplied with a titanium laboratory screw (ref. NPS VG 16 156).



Medical polymer.



Reference ETK NA 35 SB

Impression - continuation

Laboratory screws for straight, angulated and customizable abutments




m Titanium.

Version	References
Short	NPS VG 16 156
Medium	NPS VG 16 200
Long	NPS VG 16 250

Final restoration

Straight abutments



m Titanium.


Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment.

i * Do not use with Naturactis. ht = supra-implant height.

25Ncm

	ht	EP Ø 3.6	HP Ø 4.6	RP Ø 5.2	WP Ø 6
References	0.5	NPS PD 36 06	NPS PD 46 06	NPS PD 52 06	----
	1.5	NPS PD 36 16	NPS PD 46 16	NPS PD 52 16	NPS PD 60 16*
	2.5	NPS PD 36 26	NPS PD 46 26	NPS PD 52 26	NPS PD 60 26
	3.5	NPS PD 36 36	NPS PD 46 36	NPS PD 52 36	NPS PD 60 36
	4.5	NPS PD 36 46	NPS PD 46 46	NPS PD 52 46	NPS PD 60 46

7° angulated abutments



m Titanium.

Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment.

i ht = supra-implant height.

25Ncm

	ht	EP Ø 3.6	HP Ø 4.6	RP Ø 5.2	WP Ø 6
References	0.5	----	NPS PA 46 07 1	NPS PA 52 07 1	----
	1.5	NPS PA 36 07 2	NPS PA 46 07 2	NPS PA 52 07 2	NPS PA 60 07 2
	2.5	NPS PA 36 07 3	NPS PA 46 07 3	NPS PA 52 07 3	NPS PA 60 07 3
	3.5	NPS PA 36 07 4	NPS PA 46 07 4	NPS PA 52 07 4	NPS PA 60 07 4

Final restoration - continuation

15° angulated abutments




Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment. ht = supra-implant height.

Titanium.



	ht	EP Ø 3.6	NP Ø 4.6	RP Ø 5.2	WP Ø 6
References	0.5	---	NPS PA 46 15 1	NPS PA 52 15 1	---
	1.5	NPS PA 36 15 2	NPS PA 46 15 2	NPS PA 52 15 2	NPS PA 60 15 2
	2.5	NPS PA 36 15 3	NPS PA 46 15 3	NPS PA 52 15 3	NPS PA 60 15 3
	3.5	NPS PA 36 15 4	NPS PA 46 15 4	NPS PA 52 15 4	NPS PA 60 15 4

20° angulated abutments




Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment. ht = supra-implant height.

Titanium.



	ht	EP Ø 3.6	NP Ø 4.6	RP Ø 5.2	WP Ø 6
References	0.5	---	NPS PA 46 20 1	NPS PA 52 20 1	---
	1.5	NPS PA 36 20 2	NPS PA 46 20 2	NPS PA 52 20 2	NPS PA 60 20 2
	2.5	NPS PA 36 20 3	NPS PA 46 20 3	NPS PA 52 20 3	NPS PA 60 20 3
	3.5	NPS PA 36 20 4	NPS PA 46 20 4	NPS PA 52 20 4	NPS PA 60 20 4

Castable abutments with gold or chrome-cobalt base



Supplied with a titanium fixing screw (ref. NVP VPA 35) with "anti-unscrewing" treatment.



ht = supra-implant height.


Gold alloy* or chrome-cobalt alloy** + medical polymer.

ht	Gold base references	Chrome-cobalt base reference
0.5	NPS PS 46 06	NPS PCC 46 06
1.5	NPS PS 46 16	-

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 µm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

** Properties of the chrome-cobalt base and chemical composition: Chrome 26-30% / Cobalt 63-69% / Molybdenum 5-7% // Hardness (HV10): 310 // Melting range: 1370-1420°C // The temperature of the casting should not exceed 1500°C // Density: 8.3 g/cm³ // Thermal expansion: 4.1 µm / m³K // For casting, follow the manufacturer of the chrome-cobalt alloy's instructions.

Customizable abutment



Supplied with a titanium fixing screw (ref. NVP 35) with "anti-unscrewing" treatment.

Titanium.



Reference NPS PAT 46 20

CEMENTED RESTORATION ON DIRECT CLIP ABUTMENTS

For single and multi-unit prosthesis




- **Different emergence profiles and heights:**
 - Emergence diameters: 3.6 / 4.8 / 6.5 mm.
 - Coronal heights: 4 / 5.5 / 7 mm.
 - Supra-implant heights: 0.5 / 1.5 / 2.5 / 3.5 / 4.5 mm.
- **Abutments can be cut** on 2 mm at the top of the abutment (groove showing the visible limit of retouching).
- **Easy impression taking** = standardised protocol with a snap-on coping to be seated directly onto the abutment.
- **Nitrided abutments** for a better aesthetic result.
- **Laser marking** to identify emergence diameter, supra-implant height and coronal height of the abutments.
- **Kits** available include all the parts required to restore the selected abutment, the abutment should be chosen apart.
- **Colour coding of secondary components** according to the emergence profile and coronal height of the abutment:



Final abutment seating

Direct Clip abutments



*ht a = coronal height.
ht b = supra-implant height.*




Titanium.

		Ø 3.6	Ø 4.8	Ø 6.5
ht a		4	4	4
ht b	0.5	NPS PP 36 40 1	NPS PP 48 40 1	-
	1.5	NPS PP 36 40 2	NPS PP 48 40 2	NPS PP 65 40 2
	2.5	NPS PP 36 40 3	NPS PP 48 40 3	NPS PP 65 40 3
	3.5	NPS PP 36 40 4	NPS PP 48 40 4	NPS PP 65 40 4
	4.5	NPS PP 36 40 5	NPS PP 48 40 5	NPS PP 65 40 5
ht a		5.5	5.5	5.5
ht b	0.5	NPS PP 36 55 1	NPS PP 48 55 1	-
	1.5	NPS PP 36 55 2	NPS PP 48 55 2	NPS PP 65 55 2
	2.5	NPS PP 36 55 3	NPS PP 48 55 3	NPS PP 65 55 3
	3.5	NPS PP 36 55 4	NPS PP 48 55 4	NPS PP 65 55 4
	4.5	NPS PP 36 55 5	NPS PP 48 55 5	NPS PP 65 55 5
ht a		-	7	-
ht b	0.5	-	NPS PP 48 70 1	-
	1.5	-	NPS PP 48 70 2	-
	2.5	-	NPS PP 48 70 3	-
	3.5	-	NPS PP 48 70 4	-
	4.5	-	NPS PP 48 70 5	-

Temporization

Protection caps



*ht a = coronal height.
Single use.*




Medical polymer.

		Ø 3.6	Ø 4.8	Ø 6.5
ht a		4	4	4
References		APS CP 36 40	APS CP 48 40	APS CP 65 40
ht a		5.5	5.5	5.5
References		APS CP 36 55	APS CP 48 55	APS CP 65 55
ht a		-	7	-
References		-	APS CP 48 70	-

Impression

Snap-on impression copings



*Single use.
Different colors according to the version.
ht a = coronal height.*




Medical polymer.

		Ø 3.6	Ø 4.8	Ø 6.5
ht a		4	4	4
References		APS TCP 36 40	APS TCP 48 40	APS TCP 65 40
ht a		5.5	5.5	5.5
References		APS TCP 36 55	APS TCP 48 55	APS TCP 65 55
ht a		-	7	-
References		-	APS TCP 48 70	-

Impression - continuation

Snap-on open impression copings




Single use.
For modified Direct Clip abutments.



Medical polymer.

Ø	References
3.6	APS TCP 36
4.8	APS TCP 48
6.5	APS TCP 65

Direct Clip abutment analogs



Can be cut.
Different colors according to the version.
ht a = coronal height.




Titanium.

	Ø 3.6	Ø 4.8	Ø 6.5
ht a	4	4	4
References	APS H 36 40	APS H 48 40	APS H 65 40
ht a	5.5	5.5	5.5
References	APS H 36 55	APS H 48 55	APS H 65 55
ht a	-	7	-
References	-	APS H 48 70	-

Final restoration

Rotational burn-out sleeves




For multi-unit restorations.



Medical polymer, translucent white.

Ø	ht	References
3.6	10	APS BCC 36 100
4.8	10	APS BCC 48 100
6.5	7	APS BCC 65 70

Non-rotational burn-out sleeves



For single restorations.



Medical polymer, opaque white.

Ø	ht	References
3.6	10	APS BCO 36 100
4.8	10	APS BCO 48 100
6.5	7	APS BCO 65 70

Kits



Supplied without Direct Clip abutment. The kit contains a protection cap, a snap-on impression coping, a snap-on open impression coping, an analog and a burn-out sleeve corresponding to the selected abutment. ht a = coronal height.

		Ø 3.6	Ø 4.8	Ø 6.5
ht a		4	4	4
Prosthesis type	single-unit	NPS KIT E04	NPS KIT N04	NPS KIT W04
	multi-unit	NPS KIT EC4	NPS KIT NC4	NPS KIT WC4
ht a		5.5	5.5	5.5
Prosthesis type	single-unit	NPS KIT E05	NPS KIT N05	NPS KIT W05
	multi-unit	NPS KIT EC5	NPS KIT NC5	NPS KIT WC5
ht a		-	7	-
Prosthesis type	single-unit	-	NPS KIT N07	-
	multi-unit	-	NPS KIT NC7	-

SCREWED RESTORATION ON TETRA CONICAL ABUTMENTS


For single and multi-unit prosthesis



- Permits the fitting of **single prosthesis** on straight Tetra abutments only and **multi-unit prosthesis** on straight and angulated Tetra abutments, whether the implants are parallel or highly divergent.
- Design adapted to the **early or immediate loading** in the case of multiple prosthesis.
- **Straight and angulated** (17° and 30°) abutments available.
- **Wide range of supra-implant heights:** 0.5 / 1.5 / 2.5 / 3.5 / 4.5 mm.
- **Nitrided abutments** for a better aesthetic result.
- **Abutments are supplied with abutment holders** in order to facilitate their grip and placement.
- **Common secondary components for straight and angulated abutments, for multi-unit prosthesis.**
- **4.8 mm diameter shoulder** for good prosthetic support.

Final abutment seating

Straight Tetra abutments



Supplied with a single use plastic abutment holder.
For the screwing of straight abutments use the internal hexagonal key (ref. CCL HI 20 24) or the internal hexagonal mandrel (ref. UMA HI 20 26).
ht = supra-implant height.




Titanium.



35Ncm

ht	References
0.5	NPV PT 48 1
1.5	NPV PT 48 2
2.5	NPV PT 48 3
3.5	NPV PT 48 4
4.5	NPV PT 48 5

17° angulated Tetra abutments



Supplied with a titanium fixing screw (ref. NVP VPA 35) with "anti-unscrewing" treatment.



Titanium.




Supplied with a reusable titanium abutment holder. Abutment with hexagonal base (non-rotational).
ht = supra-implant height.



25Ncm

ht	References
2.5	NPV PTA 48 17 25
3.5	NPV PTA 48 17 35
4.5	NPV PTA 48 17 45

30° angulated Tetra abutments



Supplied with a titanium fixing screw (ref. NVP VPA 35) with "anti-unscrewing" treatment.



Titanium.



Supplied with a reusable titanium abutment holder. Abutment with hexagonal base (non-rotational).
ht = supra-implant height.




25Ncm

ht	References
3.5	NPV PTA 48 30 35
4.5	NPV PTA 48 30 45

Temporization

Tetra rotational temporary abutment



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.



Titanium.




For multi-unit restorations.



20Ncm

Reference UPV PMT 48 110

Non-rotational temporary abutment on Tetra straight abutment



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.



Titanium.



For single restorations.
Can be used on straight Tetra abutments only.




20Ncm

Reference UPV POT 48 110

Temporization - continuation

Tetra protection caps




Titanium.



Version	References
Short	UPV CPT 48 20
Long	UPV CPT 48 40

Impression

Tetra rotational pick-up impression coping



Supplied with a titanium laboratory screw (ref. UPV VGM 14 150).



Titanium.



For multi-unit restorations.
Open tray impression taking.



Reference UPV TM 48

Non-rotational pick-up impression coping for straight Tetra abutment



Supplied with a titanium laboratory screw (ref. UPV VGM 14 150).



Titanium.



For single restorations.
Can be used on straight Tetra abutments only.
Open tray impression taking.



Reference UPV TO 48

Pop-in impression coping




For multi-unit restorations.
Closed tray impression taking.



Titanium.

Reference UPV PI 48

Tetra abutment analog




Can be cut for use with pins.



Titanium.

Reference UPV HM 48

Laboratory guide screws





Titanium.


Version	References
Short	UPV VGM 14 40
Medium	UPV VGM 14 150
Handle + long screw	UPV VGM 14 200


Impression - continuation

Rotational Scanbody on Tetra abutment




Supplied with a titanium fixing screw (ref. NPV VG 14 105). 


 For multi-unit restorations.


 Medical polymer.


Reference ETK UN SBP

Non-rotational Scanbody on straight Tetra abutment



Supplied with a titanium fixing screw (ref. NPV VG 14 105). 


 For single restorations.


 Medical polymer.



Reference ETK UN SBPO

Final restoration

Rotational Tetra burn-out sleeve





Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment. 



 For multi-unit restorations. 

Reference UPV CMC 48 110

Tetra rotational gold base castable sleeve





Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment. 



 For multi-unit restorations. 

Reference UPV BST 48 110

Non-rotational castable sleeves on straight Tetra abutment with gold or chrome-cobalt base



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment. 

 For single restorations. 

	References
Gold base	UPV BOT 48 110
Chrome-cobalt base	UPV PCC 48 110

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 μm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

** Properties of the chrome-cobalt base and chemical composition: Chrome 26-30% / Cobalt 63-69% / Molybdenum 5-7% // Hardness (HV10): 310 // Melting range: 1370-1420°C // The temperature of the casting should not exceed 1500°C // Density: 8.3 g/cm³ // Thermal expansion: 4.1 μm / m³K // For casting, follow the manufacturer of the chrome-cobalt alloy's instructions.

SCREWED RESTORATION ON PLURAL CONICAL ABUTMENTS

For multi-unit restorations


- **For multi-unit prosthesis** on parallel or divergent implants.
- **Abutments with a 3.8 mm diameter** for thin ridges.
- **Conical bar support.**
- **Straight or angulated** (17° and 30°).
- **Wide range of supra-implant heights:** 0.5 / 1.5 / 2.5 / 3.5 / 4.5 mm.
- **Nitrided abutments** for a better aesthetic result.
- **Angulated abutments are supplied with abutment holders** in order to facilitate their picking-up and placement.



RESTORATION ON STRAIGHT PLURAL CONICAL ABUTMENTS

Final abutments seating

Plural straight abutments



ht = supra-implant height.




Titanium.



ht	References
0.5	NPV PP 31 07
1.5	NPV PP 31 18
2.5	NPV PP 31 28
3.5	NPV PP 31 38
4.5	NPV PP 31 48

Temporization

Temporary cap




Titanium.



Reference NPV CPP 38 32

Temporary abutment



Supplied with a titanium fixing screw (ref. NPV VP 18 56) with "anti-unscrewing" treatment.



Titanium.

Reference NPV PPT 46 100

Impression

Pick-up impression coping



Supplied with a titanium laboratory screw (ref. NVP VG 18 105).



Titanium.




Open tray impression taking.



Reference NPV TP 4X

Abutment analog



Can be cut for use with pins.



Titanium.

Reference NPV HP 4X

Impression - continuation

Laboratory guide screws



Titanium.

Reference NPV VG 18 105

Scanbody on Plural straight abutment



Supplied with a titanium fixing screw (ref. NPV VP 18 56) with "anti-unscrewing" treatment.




Medical polymer.

Reference ETK NA SPD

Final restoration

Burn-out sleeve



Supplied with a titanium fixing screw (ref. NPV VP 18 56) with "anti-unscrewing" treatment.




Medical polymer.

Reference NPV CCP 35 46

RESTORATION ON ANGULATED PLURAL CONICAL ABUTMENTS

Final abutments seating

17° angulated Plural abutments



Supplied with a titanium fixing screw (ref. NVP VPA 35) with "anti-unscrewing" treatment.




Supplied with a single use plastic abutment holder.
ht = supra-implant height.



Titanium.

ht	References
1.5	NPV PPA 35 17 20
2.5	NPV PPA 35 17 30
3.5	NPV PPA 35 17 40

30° angulated Plural abutments



Supplied with a titanium fixing screw (ref. NVP VPA 35) with "anti-unscrewing" treatment.



Supplied with a single use plastic abutment holder.
ht = supra-implant height.



Titanium.

ht	References
2.5	NPV PPA 35 30 30
3.5	NPV PPA 35 30 40

Temporization

Temporary cap




Titanium.



Reference NPV CPA 38 24

Temporary abutment



Supplied with a titanium fixing screw (ref. NPV VPA 14 40) with "anti-unscrewing" treatment.



Titanium.

Reference NPV TPPA 38

Impression

Pick-up impression coping



Supplied with a titanium laboratory screw (ref. NPV VPA 14 150).



Titanium.




Open tray impression taking.



Reference NPV TPA 38

Impression - continuation

Abutment analog



Can be cut for use with pins.



Titanium.

Reference NPV HPA 38


Laboratory guide screws



Titanium.

Version	References
Short	NPV VPA 14 40
Long	NPV VPA 14 150

Scanbody on Plural angulated abutment



Supplied with a titanium fixing screw (ref. NPV VPA 14 40) with "anti-unscrewing" treatment.




Medical polymer.

Reference ETK NA SPA

Final restoration

Burn-out sleeve




Supplied with a titanium fixing screw (ref. NPV VPA 14 40) with "anti-unscrewing" treatment.



Medical polymer.

Reference NPV CCPA 38

Gold base castable sleeve



Supplied with a titanium fixing screw (ref. NPV VPA 14 40) with "anti-unscrewing" treatment.



Gold alloy* + Medical polymer.

Reference NPV BSPA 38

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 µm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

REMOVABLE RESTORATION ON O-RING ABUTMENTS


For the stabilisation of removable prosthesis



- **Ideal for:**
 - total restorations,
 - restorations added to an attachment,
 - stabilisation of total restorations,
 - relining of dental prosthesis.
- Can only be used in cases where implant axes have a **maximum divergence of 15°**.
- **3 supra-implant heights:** 1.5 / 3.5 / 5.5 mm.
- **Universal ball diameter:** Ø 2.25 mm.
- **Emergence diameter:** Ø 2.9 mm.
- **3 joints with different shore hardness** available to adapt to the retention strength: 50, 60 and 70.

Final abutments seating

O-Ring abutments



For the screwing of O-Ring abutments, use the internal hexagonal key (ref. CCL HI 25 18) or the internal hexagonal mandrel (ref. CMO HI 25 26).
ht = supra-implant height.




Titanium.

ht	References
1.5	NPA OR 35 29 MALE
3.5	NPA OR 35 40 MALE
5.5	NPA OR 35 60 MALE

Impression

Impressions are taken directly on O-Ring abutments.

O-Ring abutment analog



Can be cut for use with pins.




Titanium.

Reference OPS HOBI

Final restoration

O-Ring



Supplied with an O-Ring seal of 60 shores.



Titanium + medical silicone.

Reference UPA FOR 52

O-Ring seals



Medical silicone.

Hardness	Color	References
Flexible 50 shores	Black	UPA JOR 50
Medium 60 shores	Red	UPA JOR 60
Hard 70 shores	Black	UPA JOR 70



3

etk prosthesis


Standard prosthetic components


.....

Aesthetica+²



PROSTHETIC RECOMMENDATIONS





(1) Within the limit of prosthetic axes divergences and of a reduced number of elements (8 maximum).

CEMENTED RESTORATION ON STRAIGHT AND ANGULATED ABUTMENTS


For single and multi-unit prosthesis

- **Straight or angulated** (15° and 20°) abutments available.
- **Abutments delivered with a titanium fixing screw:** the screw is secured thanks to an interior threading of the prosthetic part to avoid sudden collapse.
- Fixing screw with **“anti-unscrewing” treatment**.



Temporization and impression

Pop-in impression copings



Supplied with a titanium impression coping screw (ref. APS VG 20 140).




Closed tray impression taking.
For use as a temporary abutment, order the screw ref. APS VF 20 93.



Titanium.

Platforms	Ø a	References
4.2 (NP)	4.2	API 42 85
4.8 (RP)	4.8	API 48 85
6.5 (WP)	6.5	API 65 85

Pop-up impression copings



Supplied with a titanium impression coping screw (ref. APS VTB 20 174).



Titanium + medical polymer.



Closed tray impression taking.
Supplied with an impression cap for pick-up technique:




- Ø 4.2 : ref APS TCP 48 55
- Ø 4.8 : ref APS TCP 48 40
- Ø 6.5 : ref APS TCP 65 40

Pop-up: different colors according to the version.

Platforms	Ø a	References
4.2 (NP)	4.2	APU T 42
4.8 (RP)	4.8	APU T 48
6.5 (WP)	6.5	APU T 65

Pick-up impression copings



Supplied with a titanium impression coping screw (ref. APV VT 20 154).




Open tray impression taking.
Different colors according to the version.



Titanium.

Platforms	Ø a	References
4.2 (NP)	4.2	APE T 42
4.8 (RP)	4.8	APE T 48
6.5 (WP)	6.5	APE T 65

Implant analogs



Different colors according to the version.



Titanium.

Platforms	Ø a	References
4.2 (NP)	4.2	ALA H 42
4.8 (RP)	4.8	ALA H 48
6.5 (WP)	6.5	ALA H 65

Laboratory guide screw




Titanium.

Reference APS VG 20 140

Temporization and impression - continuation

Scanbody direct on implant



Supplied with a titanium laboratory screw (ref. APS VG 20 140).




Medical polymer.



Platforms	Ø a	References
4.2 (NP)	4.2	ETK AE NPSB
4.8 (RP)	4.8	ETK AE RPSB
6.5 (WP)	6.5	ETK AE WPSB

Bite registration abutment




Medical polymer.

Reference ETK AIEP PO

Final restoration

Straight abutments




Supplied with a titanium fixing screw (ref. APS VF 20 112) with "anti-unscrewing" treatment.



Titanium.

Platforms	References
4.2 (NP)	APS PD 48 55
4.8 (RP)	
6.5 (WP)	APS PD 65 55

15° angulated abutments




Supplied with a titanium fixing screw (ref. APS VF 20 70) with "anti-unscrewing" treatment.



Titanium.

Platforms	References
4.2 (NP)	APS PA 48 15
4.8 (RP)	
6.5 (WP)	APS PA 65 15

20° angulated abutments



Supplied with a titanium fixing screw (ref. APS VF 20 70) with "anti-unscrewing" treatment.




Titanium.

Platforms	References
4.2 (NP)	APS PA 48 20
4.8 (RP)	
6.5 (WP)	APS PA 65 20

Final restoration - continuation

Customizable abutments




Supplied with a titanium fixing screw (ref. APS VF 20 112) with "anti-unscrewing" treatment.

Titanium.

35Ncm

Platforms	Ø a	References
4.2 (NP)	4.2	APS PR 42 100
4.8 (RP)	4.8	APS PR 48 150
6.5 (WP)	6.5	APS PR 65 120

Gold base castable abutments



Supplied with a titanium fixing screw (ref. APS VF 20 70) with "anti-unscrewing" treatment.


Gold alloy* + medical polymer.

35Ncm

Platforms	References
4.2 (NP)	APS PS 48 100
4.8 (RP)	
6.5 (WP)	APS PS 65 100

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 µm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

Rotational burn-out sleeves for straight abutments




For multi-unit restorations.

Medical polymer, translucent white.

Platforms	Ø a	References
4.2 (NP)	4.2	APS CCC 42 55
4.8 (RP)	4.8	APS CCC 48 55
6.5 (WP)	6.5	APS CCC 65 55

Non-rotational burn-out sleeves for straight abutments




For single restorations.

Medical polymer, opaque white.

Platforms	Ø a	References
4.2 (NP)	4.2	APS CCO 42 55
4.8 (RP)	4.8	APS CCO 48 55
6.5 (WP)	6.5	APS CCO 65 55

Burn-out sleeves for 15 and 20° angulated abutments



Medical polymer.

Platforms	Ø a	References
4.2 (NP)	4.2	APS CA 42 47
4.8 (RP)	4.8	APS CA 48 50
6.5 (WP)	6.5	APS CA 65 70

CEMENTED RESTORATION ON DIRECT CLIP ABUTMENTS

For single and multi-unit prosthesis

- **Different coronal heights:** 4 / 5.5 / 7 mm.
- **Abutments can be retouched** on 2 mm at the top of the abutment (groove showing the visible limit of retouching).
- **Easy impression taking** = standardised protocol with a snap-on coping to be seated directly onto the abutment.
- **Kits** available include all the parts required to restore the selected abutment, the abutment should be chosen apart.
- **Colour coding of secondary components** according to the emergence profile and coronal height of the abutment:



Final abutment seating

Direct Clip abutments



Different colors according to the version.
ht = coronal height.




Titanium.



Platforms	ht	References
4.8 (RP)	4	APS PP 48 40
	5.5	APS PP 48 55
	7	APS PP 48 70
6.5 (WP)	4	APS PP 65 40
	5.5	APS PP 65 55

Temporization

Protection caps




Medical polymer, translucent white.

Platforms	ht	References
4.8 (RP)	5.8	APS CP 48 40
	7.2	APS CP 48 55
	8.8	APS CP 48 70
6.5 (WP)	5.8	APS CP 65 40
	7.2	APS CP 65 55

Impression

Snap-on impression copings




Closed tray impression taking.
Single use.
Different colors according to the version.
ht = coronal height.



Medical polymer.

Platforms	ht	References
4.8 (RP)	4	APS TCP 48 40
	5.5	APS TCP 48 55
	7	APS TCP 48 70
6.5 (WP)	4	APS TCP 65 40
	5.5	APS TCP 65 55

Snap-on open impression copings




Closed tray impression taking.
Single use.
For modified Direct Clip abutments.



Medical polymer.

Platforms	Ø a	References
4.8 (RP)	4.8	APS TCP 48
6.5 (WP)	6.5	APS TCP 65

Direct Clip abutment analogs



Can be cut for use with pins.
Different colors according to the version.
ht = coronal height.



Titanium.

Platforms	ht	References
4.8 (RP)	4	APS H 48 40
	5.5	APS H 48 55
	7	APS H 48 70
6.5 (WP)	4	APS H 65 40
	5.5	APS H 65 55

Final restoration

Rotational burn-out sleeves



For multi-unit restorations.



Medical polymer, translucent white.

Platforms	Ø a	ht	References
4.8 (RP)	4.8	10	APS BCC 48 100
6.5 (WP)	6.5	7	APS BCC 65 70

Non-rotational burn-out sleeves



For single restorations.



Medical polymer, opaque white.

Platforms	Ø a	ht	References
4.8 (RP)	4.8	10	APS BCO 48 100
6.5 (WP)	6.5	7	APS BCO 65 70

Kits



Supplied without Direct Clip abutment.

The kit contains a protection cap, a snap-on impression coping, a snap-on open impression coping, an analog, and a burn-out sleeve corresponding to the selected abutment. ht = coronal height.

		Platform	
		4.8 (RP)	6.5 (WP)
ht a		4	4
Prosthesis type	single-unit	NPS KIT N04	NPS KIT W04
	multi-unit	NPS KIT NC4	NPS KIT WC4
ht a		5.5	5.5
Prosthesis type	single-unit	NPS KIT N05	NPS KIT W05
	multi-unit	NPS KIT NC5	NPS KIT WC5
ht a		7	-
Prosthesis type	single-unit	NPS KIT N07	-
	multi-unit	NPS KIT NC7	-

SCREWED RESTORATION ON CONOCTA ABUTMENTS




∞

Aesthetica+2

For fixed screw retained prosthesis or
removable restorations on bars

PROSTHESIS - Prosthetic ranges

- The final M2 prosthetic screw, with anti-unscrewing treatment, enables **prosthesis screwing at 35 Ncm**.
- **Colour coding** of abutments and secondary components according to the emergence profile of the abutment:

<i>NP</i>	<i>RP</i>	<i>WP</i>
Ø 4.2	Ø 4.8	Ø 6.5
		



Final abutments seating

Conocta abutments



Different colors according to the version.




Titanium.



Platforms	ht	References
4.2 - 4.8 (NP - RP)	1.5	APV PI 48 15
6.5 (WP)	1.8	APV PI 65 15

Temporization

Conocta temporary abutments



Supplied with a titanium fixing screw (ref. APV VF 20 56) with "anti-unscrewing" treatment.



Titanium.



Platforms	Ø a	References
4.2 (NP)	4.2	APV PPT 42 100
4.8 (RP)	4.8	APV PPT 48 100
6.5 (WP)	6.5	APV PPT 65 100

Conocta protection caps



Supplied with a titanium fixing screw (ref. APV VF 20 56) with "anti-unscrewing" treatment.



Titanium.



Platforms	Ø a	References
4.2 (NP)	4.2	APV CP 42
4.8 (RP)	4.8	APV CP 48
6.5 (WP)	6.5	APV CP 65

Impression

Conocta pick-up impression copings



Supplied with a titanium fixing screw (ref. APE VTC 20 109).



Open tray impression taking.
Different colors according to the version.




Titanium.

Platforms	Ø a	References
4.2 (NP)	4.2	APE TC 42
4.8 (RP)	4.8	APE TC 48
6.5 (WP)	6.5	APE TC 65

Impression - continuation

Conocta abutment analogs



Can be cut for use with pins.
Different colors according to the version.



Titanium.

Platforms	Ø a	References
4.2 (NP)	4.2	ALA HC 42
4.8 (RP)	4.8	ALA HC 48
6.5 (WP)	6.5	ALA HC 65

Laboratory guide screw



Titanium.

Reference APV VG 20 150

Scanbody on Conocta abutments



Supplied with a titanium fixing screw
(ref. APV VF 20 56) with "anti-
unscrewing" treatment.



Medical polymer.



Platforms	Ø a	References
4.2 (NP)	4.2	ETK AE NPP
4.8 (RP)	4.8	ETK AE RPP
6.5 (WP)	6.5	ETK AE WPP

Final restoration

Conocta rotational burn-out sleeves



Supplied with a titanium fixing screw
(ref. APV VF 20 56) with "anti-
unscrewing" treatment.



For multi-unit restorations.



Medical polymer, translucent white.

Platforms	Ø a	References
4.2 (NP)	4.2	APV BCC 42 100
4.8 (RP)	4.8	APV BCC 48 100
6.5 (WP)	6.5	APV BCC 65 100

Conocta non-rotational burn-out sleeves



Supplied with a titanium fixing screw
(ref. APV VF 20 56) with "anti-
unscrewing" treatment.



For crowns restorations.



Medical polymer, opaque white.

Platforms	Ø a	References
4.2 (NP)	4.2	APV BCO 42 100
4.8 (RP)	4.8	APV BCO 48 100
6.5 (WP)	6.5	APV BCO 65 100

REMOVABLE RESTORATION ON O-RING ABUTMENTS

For the stabilisation of removable prosthesis



- **Ideal for:**
 - total restorations,
 - restorations added to an attachment,
 - stabilisation of total restorations,
 - relining of dental prosthesis.
- Can only be used in cases where implant axes have a **maximum divergence of 15°**.
- **2 supra-implant heights:** 2 / 4 mm.
- **Universal ball diameter:** Ø 2.25 mm.
- **Emergence diameter:** Ø 2.9 mm.
- **3 joints with different shore hardness** available to adapt to the retention strength: 50, 60 and 70.

Final abutments seating

O-Ring abutments



For the screwing of O-Ring abutments, use the internal hexagonal key (ref. CCL HI 25 18) or the internal hexagonal mandrel (ref. CMO HI 25 26).



Titanium.




Implant neck Ø	ht	References
4.2 - 4.8 - 6.5	2	APA OR 48 40
	4	APA OR 48 60

Impression

Impressions are taken directly on O-Ring abutments.

O-Ring abutment analog



Can be cut for use with pins.




Titanium.

Reference OPS HOBI

Final restoration

O-Ring



Supplied with an O'Ring seal of 60 shores.



Titanium + medical silicone.

Reference UPA FOR 52

O-Ring seals



Medical silicone.

Hardness	Color	References
Flexible 50 shores	Black	UPA JOR 50
Medium 60 shores	Red	UPA JOR 60
Hard 70 shores	Black	UPA JOR 70



3

etk prosthesis


Standard prosthetic components




Uneva+



PROSTHETIC RECOMMENDATIONS





(1) Within the limit of prosthetic axes divergences and of a reduced number of elements (8 maximum).

CEMENTED RESTORATION ON STRAIGHT AND ANGULATED ABUTMENTS

For single and multi-unit prosthesis

- **Straight and angulated** (15° and 20°).
- **3 supra-implant heights** available: 1, 2 and 3 mm
- **3 prosthesis platforms:**




- **Nitrided abutments** for a better aesthetic result.
- Abutments **delivered with a titanium fixing screw**: the screw is secured thanks to an interior threading of the prosthetic part to avoid sudden collapse.
- Fixing screw with “**anti-unscrewing**” treatment.



Temporization

Rotational temporary abutment



Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.




For multi-unit restorations.



Titanium.

Reference UPS PTR 41 100

Non-rotational temporary abutment



Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.



For single restorations.




Titanium.

Reference UPS PT 41 100

Impression

Pop-in impression copings



Supplied with a titanium laboratory screw (ref. UPS VG 20 120).




Closed tray impression taking.



Titanium.

For implant Ø	Ø	References
3.6 - 4.1	4.5	UPS TCF 50
4.8	5	UPS TCF 60

Pop-up impression coping



Supplied with a titanium impression coping screw (ref. UPS VTB 20 174).




Closed tray impression taking. Supplied with an impression cap for pick-up technique (ref. APS TCP 48 40).



Titanium + medical polymer.

Reference UPS TCU 50

Pick-up impression coping



Supplied with a titanium laboratory screw (ref. UPS VG 20 200).



Open tray impression taking.




Titanium.

Reference UPS TCO 50

Impression - continuation

Implant analogs




Can be cut for use with pins.



Titanium.

For implants Ø	Ø	References
3.6	4.1	UPS H 37
4.1		
4.8	4.9	UPS H 47

Laboratory guide screws for straight and angulated abutments



Titanium.

Lengths	References
8	UPS VG 20 80
12	UPS VG 20 120

Laboratory guide screw for castable abutments



Titanium.

Reference UPS VG 20 200

Scanbody



Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment (see page 135).




Medical polymer.



For implants Ø	Ø a	References
3.6	4.1	ETK UN 41 SB
4.1		
4.8	4.9	ETK UN 48 SB

Final restoration

Straight abutments



Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.



Titanium.



ht = supra-implant height.


* These abutments corresponding to Direct Clip abutments (see page 135).



For implants Ø	3.6 - 4.1 and 4.8 in emergency switching	4.8
Platforms Ø	4.8	5.8
References	ht 1	UPS PP 48 55 1* UPS PD 58 10
	ht 2	UPS PP 48 55 2* UPS PD 58 20
	ht 3	UPS PP 48 55 3* UPS PD 58 30

Final restoration - continuation

15° angulated abutments




Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.



Titanium.

For implants Ø	3.6 - 4.1 and 4.8 in emergence switching		4.8
Platforms Ø	4.6	5.1	5.8
References	UPS PA 46 15	UPS PA 51 15	UPS PA 58 15

20° angulated abutments




Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.



Titanium.

For implants Ø	3.6 - 4.1 and 4.8 in emergence switching		4.8
Platforms Ø	4.6	5.1	5.8
References	UPS PA 46 20	UPS PA 51 20	UPS PA 58 20

Gold castable abutments



Supplied with a titanium fixing screw (ref. UPS VCH 20 79) with "anti-unscrewing" treatment.



Gold alloy* + medical polymer.




Rotational : for multi-unit restorations.
Anti-rotational : for single-unit restorations.



For implants Ø	Ø	Rotational	Non-rotational
3.6 - 4.1 and 4.8 in emergence switching	4.6	UPS BSR 46 10	UPS BSH 46 10
4.8	5	UPS BSR 50 10	UPS BSH 50 10

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm3 // Thermal expansion: 12.4 µm / m°K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

Castable abutments



Supplied with a titanium fixing screw (ref. UPS VCH 20 79) with "anti-unscrewing" treatment.



Medical polymer.



Rotational : for multi-unit restorations.
Anti-rotational : for single-unit restorations.




For implants Ø	Ø	Rotational	Non-rotational
3.6 - 4.1 and 4.8 in emergence switching	4.6	UPS PRC 42 10	UPS PHC 42 10
4.8	5	UPS PRC 50 10	UPS PHC 50 10

CEMENTED RESTORATION ON DIRECT CLIP ABUTMENTS


For single and multi-unit prosthesis

- **Different emergence profiles and heights** are offered:
 - 2 coronal heights: 4.8 and 5.5 mm.
 - 3 gingival heights: 1, 2 and 3 mm.
- **Abutments can be retouched** on 2 mm at the top of the abutment (groove showing the visible limit of retouching).
- **Easy impression taking** = standardised protocol with a snap-on coping to be seated directly onto the abutment.
- **Nitrided abutments** for a better aesthetic result.
- **Laser marking** to identify coronal height of the abutments.
- **Kits** available include all the parts required to restore the selected abutment, the abutment should be chosen apart.
- **Colour coding of secondary components** according to the emergence profile and coronal height of the abutment:



Final abutments seating

Ø 4.8 Direct Clip abutments



Supplied with a titanium fixing screw (ref. UPS VTH 20 79) with "anti-unscrewing" treatment.



ht a = coronal height.
ht b = supra-implant height.




Titanium.

ht a		4	5.5
ht b	1	UPS PP 48 40 1	UPS PP 48 55 1
	2	UPS PP 48 40 2	UPS PP 48 55 2
	3	UPS PP 48 40 3	UPS PP 48 55 3

Temporization

Protection caps



Single use.




Medical polymer, translucent white.

ht a	References
4	APS CP 48 40
5.5	APS CP 48 55

Impression

Snap-on impression copings




Closed tray impression taking.
Single use.
Different colors according to the version.



Medical polymer.

ht a	References
4	APS TCP 48 40
5.5	APS TCP 48 55

Snap-on open impression coping




Closed tray impression taking. Single use.
For modified Direct Clip abutments.



Medical polymer.

Reference APS TCP 48

Direct Clip abutment analogs



Can be cut for use with pins.
Different colors according to the version.
ht a = coronal height.




Titanium.

ht a	References
4	APS H 48 40
5.5	APS H 48 55

Final restoration

Rotational burn-out sleeve




For multi-unit restorations.



Medical polymer, translucent white.

Reference APS BCC 48 100

Non-rotational burn-out sleeve



For single restorations.



Medical polymer, opaque white.

Reference APS BCO 48 100

Kits



Supplied without Direct Clip abutment.
The kit contains a protection cap, a snap-on impression coping, a snap-on open impression coping, an analog, and a burn-out sleeve corresponding to the selected abutment.
ht a = coronal height.

ht a	Prosthesis type	References
4	single-unit	NPS KIT N04
5.5		NPS KIT N05
4	multi-unit	NPS KIT NC4
5.5		NPS KIT NC5

SCREWED RESTORATION ON TETRA CONICAL ABUTMENTS

ω

Uneva+

PROSTHESIS - Prosthetic ranges

For single and multi-unit prosthesis

- Permits the fitting of **multi-unit prosthesis** on parallel or highly divergent implants.
- Design adapted **to the early or immediate loading** in the case of multiple prosthesis.
- **Straight and angulated** (17° and 30°) abutments available.
- **Wide range of supra-implant heights.**
- **Nitrided abutments** for a better aesthetic result.
- **Abutments are supplied with abutment holders** in order to facilitate their grip and placement.
- **Common secondary components for straight and angulated abutments, for multi-unit prosthesis.**
- **4.8 mm diameter shoulder** for good prosthetic support.



Final abutments seating

Straight Tetra abutments



Supplied with a titanium fixing screw with "anti-unscrewing" treatment.



Supplied with a single use plastic abutment holder.


For the screwing of straight abutments use the internal hexagonal key (ref. CCL HI 20 24) or the internal hexagonal mandrel (ref. UMA HI 20 26).



Titanium.

ht	Abutments references	Screws references
1	UPV PMD 41 10	UPV VMD 20 74
2	UPV PMD 41 20	UPV VMD 20 84
3	UPV PMD 41 30	UPV VMD 20 94
4	UPV PMD 41 40	UPV VMD 20 104

17° angulated Tetra abutments



Supplied with a titanium fixing screw (ref. UPV VMA 20 61) with "anti-unscrewing" treatment.




Supplied with a reusable abutment holder.



Titanium.

ht	References
2	UPV PMA 41 17 20
3	UPV PMA 41 17 30
4	UPV PMA 41 17 40

30° angulated Tetra abutments



Supplied with a titanium fixing screw (ref. UPV VMA 20 61) with "anti-unscrewing" treatment.



Supplied with a reusable abutment holder.




Titanium.

ht	References
3	UPV PMA 41 30 30
4	UPV PMA 41 30 40

Temporization

Tetra rotational temporary abutment



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.


For multi-unit restorations.

Titanium.

20Ncm

Reference UPV PMT 48 110

Non-rotational temporary abutment on Tetra straight abutment



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.

For single restorations. Can be used on straight Tetra abutments only.

Titanium.

20Ncm

Reference UPV POT 48 110

Tetra protection caps




Titanium.

10Ncm

Version	References
Short	UPV CPT 48 20
Long	UPV CPT 48 40

Impression

Tetra rotational pick-up impression coping



Supplied with a titanium laboratory screw (ref. UPV VGM 14 150).


For multi-unit restorations. Open tray impression taking.

Titanium.

5Ncm

Reference UPV TM 48

Non-rotational pick-up impression coping for straight Tetra abutment



Supplied with a titanium laboratory screw (ref. UPV VGM 14 150).


For single restorations. Can be used on straight Tetra abutments only. Open tray impression taking.

Titanium.

5Ncm

Reference UPV TO 48

Pop-in impression coping



For multi-unit restorations. Closed tray impression taking.


Titanium.

5Ncm

Reference UPV PI 48

Impression - continuation

Tetra abutment analog




Can be cut for use with pins.



Titanium.

Reference UPV HM 48

Laboratory guide screws



Titanium.

Version	References
Short	UPV VMD 14 38
Medium	UPV VGM 14 150
Handle + long screw	UPV VGM 14 200

Rotational Tetra Scanbody



Supplied with a titanium fixing screw (ref. NPV VG 14 105).




For multi-unit restorations.



Medical polymer.

Reference ETK UN SBP

Non-rotational Scanbody on straight Tetra abutment



Supplied with a titanium fixing screw (ref. NPV VG 14 105).



Medical polymer.




For single restorations.
Can be used on straight Tetra abutments only.



Reference ETK UN SBPO

Final restoration

Rotational Tetra burn-out sleeve



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.


For multi-unit restorations.

Medical polymer.

20Ncm

Reference UPV CMC 48 110

Rotational Tetra gold base castable sleeve



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.


For multi-unit restorations.

Gold alloy* + medical polymer.

20Ncm

Reference UPV BST 48 110

Non-rotational castable sleeves on straight Tetra abutment with gold or chrome-cobalt base



Supplied with a titanium fixing screw (ref. UPV VMD 14 38) with "anti-unscrewing" treatment.

For single restorations.

Gold alloy* or chrome-cobalt alloy** + medical polymer.

20Ncm


	References
Gold base	UPV BOT 48 110
Chrome-cobalt base	UPV PCC 48 110

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 μm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

** Properties of the chrome-cobalt base and chemical composition: Chrome 26-30% / Cobalt 63-69% / Molybdenum 5-7% // Hardness (HV10): 310 // Melting range: 1370-1420°C // The temperature of the casting should not exceed 1500°C // Density: 8.3 g/cm³ // Thermal expansion: 4.1 μm / m³K // For casting, follow the manufacturer of the chrome-cobalt alloy's instructions.

REMOVABLE RESTORATION ON O-RING ABUTMENTS

For the stabilisation of removable prosthesis



- **Ideal for:**
 - total restorations,
 - restorations added to an attachment,
 - stabilisation of total restorations,
 - relining of dental prosthesis.
- Can only be used in cases where implant axes have a **maximum divergence of 15°**.
- **3 supra-implant heights:** 2 / 4 / 6 mm.
- **Universal ball diameter:** Ø 2.25 mm.
- **Emergence diameter:** Ø 2.9 mm.
- **3 joints with different shore hardness** available to adapt to the retention strength: 50, 60 and 70.

Final abutments seating

O-Ring abutments



For the screwing of O-Ring abutments, use the internal hexagonal key (ref. CCL HI 25 18) or the internal hexagonal mandrel (ref. CMO HI 25 26).
ht = supra-implant height.



Titanium.




ht	References
2	UPA OR 20 MALE
4	UPA OR 40 MALE
6	UPA OR 60 MALE

Impression

Impressions are taken directly on O-Ring abutments.

O-Ring abutment analog



Can be cut for use with pins.




Titanium.

Reference OPS HOB1

Final restoration

O-Ring



Supplied with an O'Ring seal of 60 shores.



Titanium + medical silicone.

Reference UPA FOR 52

O-Ring seals



Medical silicone.

Hardness	Color	References
Flexible 50 shores	Black	UPA JOR 50
Medium 60 shores	Red	UPA JOR 60
Hard 70 shores	Black	UPA JOR 70

Ø 3 PLATFORM



1 COMMON
PROSTHETIC RANGE

1 SINGLE
CONNECTION



Naturactis and Naturall+ straight abutments (Ø 3) have a connection and prosthetic range that are not compatible with Naturactis and Naturall+ Ø 3.5 - 4 - 4.5 - 5 implants.

3

etk prosthesis

Standard prosthetic components

.....

Naturactis Ø 3
Naturall+ Ø 3

PROSTHETIC RECOMMENDATIONS

Temporary restorations

Final restorations

Single

Single



Cemented



20Ncm

Temporary abutments
non-rotational
p 149



20Ncm

Straight and angulated abutments
p 150



20Ncm

Castable abutments
gold base
p 150



30Ncm

Direct Clip abutments
p 152



20Ncm

Esthetibase interfaces
non-rotational
p 163



20Ncm

Straight and angled abutments
p 150



30Ncm

O-Ring abutments
p 155



20Ncm

Castable abutments
gold base
p 150



30Ncm

Direct Clip abutments
p 152



20Ncm

Esthetibase interfaces
rotational
p 163

(1) Within the limit of prosthetic axes divergences and of a reduced number of elements (8 maximum).

CEMENTED RESTORATION ON STRAIGHT AND ANGULATED ABUTMENTS

For single and multi-unit prosthesis
(on lower incisors)

- **1 single prosthetic platform:**






- **3 supra-implant heights:** 1 / 3 / 5 mm.
- **Straight or angulated** (7° and 15°) abutments available.
- **Nitrided abutments** for a better aesthetic result.
- **Abutments delivered with a titanium fixing screw:** the screw is secured thanks to an interior threading of the prosthetic part to avoid sudden collapse.
- Fixing screw with “**anti-unscrewing**” treatment.
- **Laser marking** to identify the prosthesis platform and supra-implant height of the abutments.




Temporization

Non-rotational temporary abutment




 Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.
  20Nm



 Titanium.



Reference NPS PPT 30 1

Impression

Pop-in impression coping






 Supplied with a titanium impression coping screw (ref NPS VTB 14 156).
  Titanium.



 Closed tray impression taking.
  5 Nm

Reference NPI 31

Pick-up impression coping





 Supplied with a titanium laboratory screw (ref NPS VG 14 200).
  Titanium.


 Open tray impression taking.
  5 Nm

Reference NPET 30

Implant analog





 Can be cut for use with pins.

 Titanium.

Reference NLA H 30

Laboratory guide screws



 Titanium.

Version	References
Medium	NPS VG 14 200
Long	NPS VG 14 250

Impression - continuation

Scanbody direct on implant



Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.




Medical polymer.

Reference ETK NA 30 SB

Final restoration

Straight abutments



Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.



ht = supra-implant height.



Titanium.

ht	References
0.5	NPS PD 30 06
2.5	NPS PD 30 26
4.5	NPS PD 30 46

7° angulated abutments



Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.



ht = supra-implant height.



Titanium.

ht	References
0.5	NPS PA 30 07 1
2.5	NPS PA 30 07 3

15° angulated abutments



Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.




ht = supra-implant height.



Titanium.

ht	References
0.5	NPS PA 30 15 1
2.5	NPS PA 30 15 3

Gold base castable abutment



Supplied with a titanium fixing screw (ref. NVP 31) with "anti-unscrewing" treatment.



Gold alloy* + medical polymer.



Reference NPS PS 30 16

* Properties of the gold base and chemical composition: Gold (Au) 58.25% +/- 1% // Platinum (Pt) 21.90% +/- 1% // Palladium (Pd) 19.41% +/- 1% // Iridium (Ir) 0.44% + 0.5%/- 0% // Hardness (HV) > 160 // Solidus - Liquidus: 1400-1490°C // Density: 17.5 g/cm³ // Thermal expansion: 12.4 µm / m³K // Choose a casting alloy in line with ISO 9693, ISO 1891 and ISO 1562 standards that is compatible with a melting point under 1350°C // Source: Ceramicor® - Cendres & Métaux

CEMENTED RESTORATION ON DIRECT CLIP ABUTMENTS


For single and multi-unit prosthesis
(on lower incisors)



- **3 heights:** 1 / 3 / 5 mm.
- **Abutments can be retouched** on 2 mm at the top of the abutment (groove showing the visible limit of retouching).
- **Easy impression taking** = standardised protocol with a snap-on coping to be seated directly onto the abutment.
- **Nitrided abutments** for a better aesthetic result.
- **Laser marking** to identify emergence diameter, supra-implant height and coronal height of the abutments.
- **Kits** available include all the parts required to restore the selected abutment, the abutment should be chosen apart.

Final abutments seating

Direct Clip abutments



ht = supra-implant height.




Titanium.



Ø	ht	Coronal ht	References
3.6	1	5.5	NPS PP 30 55 1
	3		NPS PP 30 55 3
	5		NPS PP 30 55 5

Temporization

Protection cap



Single use.




Medical polymer, opaque white.

Reference APS CP 36 55

Impression

Snap-on impression coping



*Closed tray impression taking.
Single use.
ht = coronal height.*



Medical polymer.

Reference APS TCP 36 55

Snap-on open impression coping




*Closed tray impression taking.
Single use.
For modified Direct Clip abutments.*



Medical polymer.

Reference APS TCP 36

Direct Clip abutment analog



Can be cut for use with pins.




Titanium.

Reference APS H 36 55

Final restoration

Rotational burn-out sleeve




For multi-unit restorations.



Medical polymer, translucent white.

Reference APS BCC 36 100

Non-rotational burn-out sleeve



For single restorations.



Medical polymer, opaque white.

Reference APS BCO 36 100

Kits



Supplied without Direct Clip abutment. The kit contains a protection cap, a snap-on impression coping, a snap-on open impression coping, an analog, and a burn-out sleeve corresponding to the selected abutment.

Prosthesis type	References
single-unit	NPS KIT E05
multi-unit	NPS KIT EC5

REMOVABLE RESTORATION ON O-RING ABUTMENTS

For the stabilisation of removable prosthesis



- **Ideal for:**
 - total restorations,
 - restorations added to an attachment,
 - stabilisation of total restorations,
 - relining of dental prosthesis.
- Can only be used in cases where implant axes have a **maximum divergence of 15°**.
- **Universal ball diameter:** Ø 2.25 mm.
- **Emergence diameter:** Ø 2.9 mm.
- **3 joints with different shore hardness** available to adapt to the retention strength: 50, 60 and 70.

Final abutments seating

O-Ring abutments



For the screwing of O-Ring abutments, use the internal hexagonal key (ref. CCL HI 25 18) or the internal hexagonal mandrel (ref. CMO HI 25 26).



Titanium.




Reference NPA OR 31 29 MALE

Impression

Impressions are taken directly on O-Ring abutments.

O-Ring abutment analog



Can be cut for use with pins.




Titanium.

Reference OPS H0BI

Final restoration

O-Ring



Supplied with an O'Ring seal of 60 shores.



Titanium + medical silicone.



Reference UPA FOR 52

O-Ring seals



Medical silicone.

Hardness	Color	References
Flexible 50 shores	Black	UPA JOR 50
Medium 60 shores	Red	UPA JOR 60
Hard 70 shores	Black	UPA JOR 70



3

etk prosthesis

Standard prosthetic components



Obi Ø 2.7



REMOVABLE RESTORATION



For the stabilisation of removable prosthesis




- **Overdenture stabilization on thin crests.**
- **Universal ball diameter:** Ø 2.25 mm.
- **Emergence diameter:** Ø 2.9 mm.
- **3 joints with different shore hardness** available to adapt to the retention strength: 50, 60 and 70.

Impression

Impressions are taken directly on implants.

O-Ring abutment analog



2.25
15


i Can be cut for use with pins.

m Titanium.

Reference OPS H0BI

Final restoration

O-Ring




5.2

i Supplied with an O'Ring seal of 60 shores.

m Titanium + medical silicone.


Reference UPA FOR 52

O-Ring seals



m Medical silicone.

Hardness	Color	References
Flexible 50 shores	Black	UPA JOR 50
Medium 60 shores	Red	UPA JOR 60
Hard 70 shores	Black	UPA JOR 70



3

etk prosthesis

CAD-CAM customised prosthesis




- 160 Titanium Esthetibase interfaces
- 162 Scanbody
- 164 CAD-CAM works



TITANIUM ESTHETIBASE INTERFACES

For IPS e.max® or zirconia
CAD-CAM prosthesis

- **Titanium interfaces** for:
 - CAD-CAM prosthesis.
 - Profile Designer iphysio® (see pages 80-81).
- **Nitrided surface** = invisible under the gingiva.
- **2 bases:**
 - Rotational for multi-unit cases.
 - Non-rotational for single cases.
- **3 platforms :**



Cerec® compatibility

Esthetibase interfaces compatible with etk implants

On implant

<i>For suprastructure on implant</i>	<i>Platform</i>	<i>For crown</i>	<i>For bridge</i>
<i>Naturall+ Ø 3, Naturactis Ø 3</i>	S	NPC PTO 30 16	-
<i>Naturall+, Natea+, Naturactis, Ht 1</i>	S	NPC PTO 40 06	NPC PTC 40 06
<i>Naturall+, Natea+, Naturactis, Ht 2</i>	S	NPC PTO 40 16	NPC PTC 40 16
<i>Aesthetica+² NP Ø 4.2</i>	L	APC PTO 42	APC PTC 42
<i>Aesthetica+² RP Ø 4.8</i>	XL	APC PTO 48	APC PTC 48
<i>Aesthetica+² WP Ø 6.5</i>	-	APC PTO 65	APC PTC 65
<i>Uneva+ Ø 3.6 and 4.1</i>	L	UPC PTO 41	UPC PTC 41
<i>Uneva+ Ø 4.8 emergence switching Ø 4.1 if low angulation</i>			

On conical abutment

<i>For suprastructure on conical abutment</i>	<i>Platform</i>	<i>For crown</i>	<i>For bridge</i>
<i>Tetra Esthetibase</i>	-	-	UPV PTC

SCANBODY

For greater precision during scanning

Scanbodies are digital impression coping devices used to determine the position of the implant with regards to the rest of the mouth (digital intra-oral scanning) or on a master cast model (laboratory scanning system).


Their specific shape guarantees very high precision during scanning.

After scanning and recognition of the Scanbody, the design process can begin.

Scanbodies can be sterilised in an autoclave for intra-oral use.

They are compatible with the following softwares:

- Dental Wings®
- Imetric® - Exocad®
- 3Shape®



etk Scanbodies



etk scanbodies are sold by unit or by 6.

Direct on implant

<i>For suprastructure on implant</i>	<i>Laser marking</i>	<i>References</i>
<i>Naturall+ Ø 3 / Naturactis Ø 3</i>	NAT 3.0	ETK NA 30SB
<i>Naturall+, Natea+, Naturactis</i>	NAT 3.5	ETK NA 35SB
<i>Aesthetica+² NP Ø 4.2</i>	AEST NP	ETK AE NPSB
<i>Aesthetica+² RP Ø 4.8</i>	AEST RP	ETK AE RPSB
<i>Aesthetica+² WP Ø 6.5</i>	AEST WP	ETK AE WPSB
<i>Uneva+ Ø 3.6 and Ø 4.1</i>	UNE 4.1	ETK UN 41SB
<i>Uneva+ Ø 4.8</i>	UNE 4.8	ETK UN 48SB

On conical abutment

<i>For suprastructure on abutment</i>	<i>Laser marking</i>	<i>References</i>	
<i>Naturall+, Natea+, Naturactis</i>	<i>Straight Plural</i>	PLURAL DRT	ETK NA SPD
	<i>Angulated Plural</i>	PLURAL ANG	ETK NA SPA
<i>Aesthetica+² NP Ø 4.2</i>	<i>Conocta</i>	AE SBP N	ETK AE NPP
<i>Aesthetica+² RP Ø 4.8</i>		AE SBP R	ETK AE RPP
<i>Aesthetica+² WP Ø 6.5</i>		AE SBP W	ETK AE WPP
<i>Naturall+, Natea+, Naturactis Uneva+</i>	<i>Rotational Tetra</i>	TETRA	ETK UN SBP
	<i>Non-rotational Tetra (only on straight Tetra abutment)</i>	TETRA	ETK UN SBPO

On Esthetibase interface

<i>For suprastructure on implant</i>	<i>Laser marking</i>	<i>References</i>
<i>Platform S</i>	E-BASE S	ETK ESTH S SB
<i>Platform L or XL</i>	E-BASE L	ETK ESTH L SB
<i>Tetra</i>	TETRA	ETK ESTH X SB

CAD-CAM WORKS

Our expert centre, based in France, can design and manufacture single prosthesis or bridges, on natural teeth or implants, as well as customised abutments and bars using CAD and CAM.



Customised abutments

- Implant connections made on 11 axes CNN machines with a precision of 5 µm to guarantee the precision of implant/abutment assembly and sealing.
- Enables perfect defining of the cervical contour.
- Enables treatment of cases that cannot be treated using standard abutments.
- Designed to be compatible with the future tooth and sleeve.
- Enables perfect alignment of abutments.



Titanium



Zirconium or
IPS e.max® on
Esthetibase

Simple or anatomical bars and bridges on abutments or direct implants

- High precision.
- Finishing quality.
- Complies with the Sheffield test.



Zirconium bridge
on Esthetibase



Anatomical bar

For suprastructure	Customised abutments		Bars	
	on Esthetibase interfaces	titanium	on implants	on abutments
Naturall+ Ø 3, Naturactis Ø 3	x	x	-	-
Naturall+, Natea+, Naturactis	x	x	-	x
Aesthetica+² NP Ø 4.2	x	x	-	x
Aesthetica+² RP Ø 4.8	x	x	x (titanium)	x
Aesthetica+² WP Ø 6.5				
Uneva+ Ø 3.6 and Ø 4.1	x	x	x	x
Uneva+ Ø 4.8 (switching emergence Ø 4.1 if low angulation)	x	x	x	x


x : Available
- : Not available



3

etk prosthesis

All in bar® system



ALL^{IN}BAR[®] SYSTEM

Create a rigid framework for a final bridge on the same day as placing implants!


Objectives

Create an extremely rigid titanium framework/bar:


- to distribute the efforts,
- to stop all micromovements,
- to strengthen the bridge and make it "definitive".

The framework obtained has all the mechanical qualities of a machined framework at a highly competitive cost for creation in 6 hours.


1 - Placing implants and screwing Tetra abutments



2 - Occlusion registration and impression taking




3 - Screwing of final bridge



Winged copings

Possibility of retouching wings by cutting the height and/or length.

Large holes for a very high retention of the resin in the framework to prevent cosmetic fractures.




2 wings that are flexible horizontally and very rigid vertically. The horizontal flexibility allows the alloys to be curved to fit the patient's dental arch.

Mounted onto Tetra abutments of our implant ranges.


Sandblasted coping to allow the resin to have a good grip.

Shifted winged coping



For cases where the implant is not centered on the axis of the teeth.


Centered winged coping



Slight offsetting for overlapping the copings between the abutments, for cases where the implant is placed perfectly on the axis of the teeth.

Winged copings

Centered winged copings




Supplied with a fixing screw (ref. UPV VMD 14 38).

Titanium.

20Ncm

Version	References
Short	ARS PT AC
Long	ARS PT AC 30

Shifted winged copings



Supplied with a fixing screw (ref. UPV VMD 14 38).


Titanium.

20Ncm

Version	References
Short	ARS PT AT
Long	ARS PT AT 30

Impression and temporization

Short pick-up impression coping




Supplied with a fixing screw (ref. UPV VGM 14 100).

Titanium.

Batch of 4 pieces.

Reference ARS TT 74 4

Occlusion abutment




For occlusion taking by indirect method on Tetra abutment.

Medical polymer.

Quantity	References
Unité	UPV P0 48
Lot de 10	UPV P0 48 10


Tetra abutment caps



Titanium.

Version	References
Short	UPV CPT 48 20
Long	UPV CPT 48 40

Storage kit



This kit of 2 trays allows you to group together all of your useful components for All in bar® permanent resin/titanium bridges. It will accompany your cases between you and your dental technician.

Tray 2 (one space per component to help you manage resupplying):

- 8 short centered winged copings
- 2 long centered winged copings
- 4 short shifted winged copings
- 2 long shifted winged copings
- 16 winged coping screws

Tray 1:

- 8 short pick-up impression copings
- 8 occlusion abutments
- 8 short Tetra abutment caps

Reference ARSK 74



3

etk prosthesis

Common instruments




- 174 Prosthetic torque wrench
- 175 Prosthetic kit
- 176 Keys & mandrels



COMMON INSTRUMENTS

Prosthetic torque wrench



For precision tightening of prosthetic parts.
If the prosthetic rehabilitations are not assembled with the recommended tightening torque (see user guides) the screw may come loose after a few solicitations.

Furthermore,

- a tightening torque too low will not create enough strain to absorb the stress suffered by the screw.
- a tightening torque too high could damage the implant thread or deform the screw not allowing it to absorb the masticatory stress.

Warning: The torque values are indicated for permanent prosthesis. In case of immediate loading, we recommend using a lower torque value and tightening to the definitive torque later on.

- Made of surgical stainless steel.
- Detachable for easier cleaning.
- Different torque adjustments available: 10, 15, 20, 25, 30, 35 and 40 Ncm.
- Automatic stop: tightening stops automatically when the pre-set value is reached.

Store the key at a torque less than 10 Ncm to avoid premature ageing.


Reference CCC 35



Download the cleaning instructions



Option



The torque wrench can be used with all implant systems thanks to the adapter key. The adapter can be used only on prosthetic systems with a torque of less than or equal to 35 Ncm.

Reference CAD 115

Prosthetic kit



Necessary instruments for the screwing of our prosthetic parts in all implant systems.


Reference TPK 00 P6

Contents of the kit reference TPK 00 P6

Instruments	Prosthetic torque wrench 10 - 40 Ncm		CCC 35
	External hexagonal keys	<i>short</i>	CCL HE 12 18
		<i>medium</i>	CCL HE 12 22
		<i>long</i>	CCL HE 12 30
	External hexagonal mandrels	<i>short</i>	CMA HE 12 22
		<i>long</i>	CMA HE 12 26
	Internal hexagonal keys	<i>for straight Tetra abutment</i>	CCL HI 20 24
		<i>short - for Obi implant and O-Ring abutment</i>	CCL HI 25 18
		<i>long - for Obi implant and O-Ring abutment</i>	CCL HI 25 26
	Internal hexagonal mandrels	<i>for straight Tetra abutment</i>	UMA HI 20 26
<i>for Obi implant and O-Ring abutment</i>		CMO HI 25 26	

Keys and mandrels


External hexagonal keys



A key for all our prosthetic parts, except the O-Ring abutments of removable restoration and the conical Tetra abutments.

Version	References
Short	CCL HE 12 18
Medium	CCL HE 12 22
Long	CCL HE 12 30


Internal hexagonal key



For straight Tetra abutment.

Reference CCL HI 20 24


O-Ring internal hexagonal keys



For Obi implant and O-Ring abutment.


Version	References
Short	CCL HI 25 18
Long	CCL HI 25 26

External hexagonal mandrels



Version	References
Short	CMA HE 12 22
Long	CMA HE 12 26

Internal hexagonal mandrel



For straight Tetra abutment.

Reference UMA HI 20 26


Keys and mandrels - continuation

O-Ring internal hexagonal mandrel



For Obi implant and O-Ring abutment.

Reference CMO HI 25 26



Part 4

Other CAD-CAM prosthesis brands



- 180 Titanium Esthetibase interfaces
- 182 Scanbody
- 184 CAD-CAM works

TITANIUM ESTHETIBASE INTERFACES


For IPS e.max® or zirconia CAD-CAM prosthesis



Screw heads compatible with etk external hexagonal keys:

- short: *reference* CCL HE 12 18
- medium: *reference* CCL HE 12 22
- long: *reference* CCL HE 12 30

- **Titanium interfaces for:**
 - CAD-CAM prosthesis.
 - Profile Designer iphysio® (see pages 80-81).
- **Nitrided surface** = invisible under the gingiva.
- **2 bases:**
 - Rotational for multi-unit cases.
 - Non-rotational for single cases.
- **3 platforms :**



Cerec® compatibility

Other brands compatible Esthetibase interfaces

Brand	For suprastructure on implant	Platform	For crown	For bridge
AB DENTAL®	implant I® platform Ø 3.75 compatible	L	ZIM_SC.35.PTO	ZIM_SC.35.PTC
ANTHOGYR®	Axiom® compatible	L	ANT_PTO.AX	-
	Anthofit® HE Ø 4.1 compatible	L	UPC_PTO.41	UPC_PTC.41
	Anthofit® HE Ø 5 compatible	XL	BIO_EX.50.PTO	BIO_EX.50.PTC
	Ossfit® platform Ø 4.8 compatible	XL	APC_PTO.48	APC_PTC.48
	Ossfit® platform Ø 6.5 compatible	/	APC_PTO.65	APC_PTC.65
ASTRA®	Ø 3.5 and 4 connection Ocean® Ht1 compatible	S	NPC_PTO.40.06	NPC_PTC.40.06
	Ø 3.5 and 4 connection Ocean® Ht2 compatible	S	NPC_PTO.40.16	NPC_PTC.40.16
	Ø 4.5 and 5 connection Lilas® compatible	XL	NPC_PTO.50	NPC_PTC.50
BIOHORIZONS®	yellow connection Ø 3.5 compatible	L	ZIM_SC.35.PTO	ZIM_SC.35.PTC
	connection green Ø 4.5 compatible	L	ZIM_SC.45.PTO	ZIM_SC.45.PTC
	connection blue Ø 5.7 compatible	/	ZIM_SC.57.PTO	ZIM_SC.57.PTC

Other brands compatible Esthetibase interfaces - continuation

<i>Brand</i>	<i>For suprastructure on implant</i>	<i>Platform</i>	<i>For crown</i>	<i>For bridge</i>
BIOMET 3i®	<i>Certain® Ø 3.4 purple compatible</i>	S	BIO_CE.3.PTO	BIO_CE.3.PTC
	<i>Certain® Ø 4.1 blue compatible</i>	S	BIO_CE.4.PTO	BIO_CE.4.PTC
	<i>Certain® Ø 5 yellow compatible</i>	XL	BIO_CE.5.PTO	BIO_CE.5.PTC
	<i>Certain® Ø 6 green compatible</i>	/	BIO_CE.6.PTO	BIO_CE.6.PTC
	<i>External Hex. Ø 4.1 compatible</i>	L	UPC_PTO.41	UPC_PTC.41
	<i>External Hex. Ø 5 compatible</i>	XL	BIO_EX.50.PTO	BIO_EX.50.PTC
	<i>External Hex. Ø 6 compatible switching emergence Ø 5 if low angulation</i>	XL	BIO_EX.50.PTO	BIO_EX.50.PTC
BIOTECH®	<i>Kontakt® (Ø 3.6 - 4.2 - 4.8 - 5.4) compatible</i>	S	BTH_KO.X.PTO	-
EASY IMPLANT®	<i>connection Ocean® Ht1 compatible</i>	S	NPC_PTO.40 06	NPC_PTC.40.06
	<i>connection Ocean® Ht2 compatible</i>	S	NPC_PTO.40 16	NPC_PTC.40.16
	<i>connection Lilas® compatible (emergence Ø 5.1) Ht1</i>	L	NPC_PTO.4X	NPC_PTC.4X
	<i>connection Lilas® compatible (emergence Ø 5.8) Ht1</i>	XL	NPC_PTO.4X.58	NPC_PTC.4X.58
	<i>connection Lilas® compatible (emergence Ø 5.8) Ht2</i>	XL	NPC_PTO.50	NPC_PTC.50
GLOBAL D®	<i>Global D® Tekka® Inkone® compatible</i>	L	TEK_PTO.INK	-
	<i>Global D® Tekka® Hexa Color® green compatible</i>	S	TEK_PTO.HE.VE	-
	<i>Global D® Tekka® Hexa Color® orange compatible</i>	S	TEK_PTO.HE.OR	-
	<i>Global D® Tekka® Hexa Color® yellow compatible</i>	L	TEK_PTO.HE.JA	TEK_PTC.HE.JA
	<i>Global D® Tekka® Hexa Color® purple compatible</i>	XL	TEK_PTO.HE.VI	TEK_PTC.HE.VI
	<i>Global D® Serf® EVL® Ø 3.3 yellow compatible</i>	S	SER_EV.3.PTO	SER_EV.3.PTC
	<i>Global D® Serf® EVL® Ø 4 blue compatible</i>	L	SER_EV.4.PTO	SER_EV.4.PTC
IMPLANT DIRECT®	<i>Legacy® green compatible</i>	L	ZIM_SC.35.PTO	ZIM_SC.35.PTC
	<i>Legacy® purple compatible</i>	L	ZIM_SC.45.PTO	ZIM_SC.45.PTC
	<i>Legacy® yellow compatible</i>	/	ZIM_SC.57.PTO	ZIM_SC.57.PTC
	<i>Interactive® purple compatible</i>	S	NOB_AC.N.PTO	NOB_AC.N.PTC
	<i>Interactive® yellow compatible</i>	L	NOB_AC.R.PTO	NOB_AC.R.PTC
MIS®	<i>Seven® M4 SP purple compatible</i>	L	ZIM_SC.35.PTO	ZIM_SC.35.PTC
	<i>Seven® M4 WP green compatible</i>	L	ZIM_SC.45.PTO	ZIM_SC.45.PTC
NOBEL BIO CARE®	<i>Active® Ø 3.5 NP et Replace® CC compatible</i>	S	NOB_AC.N.PTO	NOB_AC.N.PTC
	<i>Active® Ø 4.3 and 5 RP - Replace® CC compatible</i>	L	NOB_AC.R.PTO	NOB_AC.R.PTC
	<i>Replace® NP red compatible</i>	L	NOB_RE.N.PTO	NOB_RE.N.PTC
	<i>Replace® RP yellow compatible</i>	L	NOB_RE.R.PTO	NOB_RE.R.PTC
	<i>Replace® WP blue compatible</i>	XL	NOB_RE.W.PTO	NOB_RE.W.PTC
	<i>Branemark® NP compatible</i>	S	NOB_BR.N.PTO	NOB_BR.N.PTC
	<i>Branemark® RP compatible</i>	L	UPC_PTO.41	UPC_PTC.41
STRAUMANN®	<i>BoneLevel® NC Ø 3.3 compatible</i>	S	STR_BL.N.PTO	STR_BL.N.PTC
	<i>BoneLevel® RC Ø 4.1 and 4.8 compatible</i>	S	STR_BL.R.PTO	STR_BL.R.PTC
	<i>Tissue level® RN compatible</i>	XL	APC_PTO.48	APC_PTC.48
	<i>Tissue level® WN compatible</i>	/	APC_PTO.65	APC_PTC.65
ZIMMER®	<i>Tapered Screw-Vent® Ø 3.5 compatible</i>	L	ZIM_SC.35.PTO	ZIM_SC.35.PTC
	<i>Tapered Screw-Vent® Ø 4.5 compatible</i>	L	ZIM_SC.45.PTO	ZIM_SC.45.PTC
	<i>Tapered Screw-Vent® Ø 5.7 compatible</i>	/	ZIM_SC.57.PTO	ZIM_SC.57.PTC

SCANBODY

For greater precision during scanning



Screw heads compatible with etk external hexagonal keys:

- short: *reference* CCL HE 12 18
- medium: *reference* CCL HE 12 22
- long: *reference* CCL HE 12 30

Scanbodies are digital impression coping devices used to determine the position of the implant with regards to the rest of the mouth (digital intra-oral scanning) or on a master cast model (laboratory scanning system).

Their specific shape guarantees very high precision during scanning.

After scanning and recognition of the Scanbody, the design process can begin.

Scanbodies can be sterilised in an autoclave for intra-oral use.

They are compatible with the following softwares:

- Dental Wings®
- Imetric® - Exocad®
- 3Shape®

Other brands compatible Scanbodies



etk scanbodies are sold by unit or by 6.

Brand	Connection type	For suprastructure	Laser marking on Scanbody	Etk references
AB DENTAL®	Direct implant	implant I® platform Ø 3.75 compatible	ZIM SC 35	ZIM_SC.35.SB
ANTHOGYR®	Direct implant	Axiom® compatible	AXIOM	ANT_AX.SB
		Anthofit® HE Ø 4.1 compatible	UNE 4,1	ETK_UN.41.SB
		Anthofit® HE Ø 5 compatible	UNE 4,1	ETK_UN.41.SB
		Ossfit® Ø 4.8 compatible	AEST RP	ETK_AE.RPSB
		Ossfit® Ø 6.5 compatible	AEST WP	ETK_AE.WPSB

Other brands compatible Scanbodies - continuation

<i>Fabricant</i>	<i>Connection type</i>	<i>For suprastructure</i>	<i>Laser marking on Scanbody</i>	<i>Etk references</i>
ASTRA®	<i>Direct implant</i>	<i>Ø 3.5 and Ø 4 connection Ocean compatible</i>	NAT 3,5	ETK_NA.35SB
		<i>Ø 4.5 and Ø 5 connection Lilas compatible</i>	NAT 5,0	ETK_NA.50SB
BIOHORIZONS®	<i>Direct implant</i>	<i>connection yellow Ø 3.5 compatible</i>	ZIM SC 35	ZIM_SC.35.SB
		<i>connection green Ø 4.5 compatible</i>	ZIM SC 45	ZIM_SC.45.SB
		<i>connection blue Ø 5.7 compatible</i>	ZIM SC 57	ZIM_SC.57.SB
BIOMET 3i®	<i>Direct implant</i>	<i>Certain Ø 3.4 - purple</i>	BIO 3i CER 3,4	BIO_CER.34.SB
		<i>Certain Ø 4.1 - Ø 4.1 blue / Ø 5 yellow / Ø 6 green compatible</i>	BIO 3i CER	BIO_CER.SB
		<i>External Hex. Ø 4 compatible</i>	UNE 4,1	ETK_UN.41SB
		<i>External Hex. Ø 5 compatible</i>	UNE 4,1	ETK_UN.41SB
	<i>External Hex. Ø 6 compatible</i>	UNE 4,1	ETK_UN.41SB	
	<i>On conical abutment</i>	<i>Low Profil - Ø 4.8 compatible</i>	TETRA	ETK_UN.SBP
BIOTECH®	<i>Direct implant</i>	<i>Kontakt Ø 3.6 - Ø 4.2 - Ø 4.8 - Ø 5.4 compatible</i>	BIOTECH K	BTH_KO.X.SB
EASY IMPLANT®	<i>Direct implant</i>	<i>narrow platform Ocean compatible</i>	NAT 3,5	ETK_NA.35SB
		<i>large platform Lilas compatible</i>	NAT 4X	ETK_NA.4XSB
GLOBAL D® (TEKKA® & SERF®)	<i>Direct implant</i>	<i>EVL Ø 3.3 yellow compatible</i>	SERF EV33	SER_EV.3.SB
		<i>EVL Ø 4 blue - Ø 5 pink compatible</i>	SERF EV 4	SER_EV.4.SB
		<i>Inkone® compatible</i>	INKONE	TEK_SB.INK
		<i>Hexa Color® green - orange compatible</i>	HEXA OR	TEK_SB.HEX.OR
		<i>Hexa Color® yellow - purple - blue - grey compatible</i>	HEXA JA	TEK_SB.HEX.JA
IMPLANT DIRECT®	<i>Direct implant</i>	<i>Legacy® green compatible</i>	ZIM SC 35	ZIM_SC.35.SB
		<i>Legacy® purple compatible</i>	ZIM SC 45	ZIM_SC.45.SB
		<i>Legacy® yellow compatible</i>	ZIM SC 57	ZIM_SC.57.SB
		<i>Interactive® purple compatible</i>	ACTIV NP	NOB_AC.N.SB
		<i>Interactive® yellow compatible</i>	ACTIV RP	NOB_AC.N.SB
MIS®	<i>Direct implant</i>	<i>Seven® M4 SP purple compatible</i>	ZIM SC 35	ZIM_SC.35.SB
		<i>Seven® M4 WP green compatible</i>	ZIM SC 45	ZIM_SC.45.SB
NOBEL®	<i>Direct implant</i>	<i>Active® Ø 3.5 NP and Replace® CC compatible</i>	ACTIV NP	NOB_AC.N.SB
		<i>Active® Ø 4.3 and Ø 5 RP - Replace® CC compatible</i>	ACTIV RP	NOB_AC.R.SB
		<i>Replace® NP red compatible</i>	REPLACE NP	NOB_RE.N.SB
		<i>Replace® RP yellow compatible</i>	REPLACE RP	NOB_RE.R.SB
		<i>Replace® WP blue compatible</i>	REPLACE WP	NOB_RE.W.SB
		<i>Brånemark® NP compatible</i>	BRAN NP	NOB_BR.N.SB
		<i>Brånemark® RP compatible</i>	UNE 4,1	ETK_UN.41SB
	<i>Brånemark® WP compatible</i>	BRAN NP	NOB_BR.W.SB	
	<i>On conical abutment</i>	<i>compatible Multi-unit® RP compatible</i>	TETRA	ETK_UN.SBP
STRAUMANN®	<i>Direct implant</i>	<i>Tissue Level RN® compatible</i>	AEST RP	ETK_AE.RPSB
		<i>Tissue Level WN® compatible</i>	AEST WP	ETK_AE.WPSB
		<i>Bone Level NC® compatible</i>	Bone NC	STR_BL.N.SB
		<i>Bone Level RC® compatible</i>	Bone RC	STR_BL.R.SB
ZIMMER®	<i>Direct implant</i>	<i>Tapered Screw-Vent® Ø 3.5 green compatible</i>	ZIM SC 35	ZIM_SC.35.SB
		<i>Tapered Screw-Vent® Ø 4.5 purple compatible</i>	ZIM SC 45	ZIM_SC.45.SB
		<i>Tapered Screw-Vent® Ø 5.7 yellow compatible</i>	ZIM SC 57	ZIM_SC.57.SB

CAD-CAM WORKS

Our expert centre, based in France, can design and manufacture single prosthesis or bridges, on natural teeth or implants, as well as customised abutments and bars using CAD and CAM.

Screw heads compatible with etk external hexagonal keys:

- short: *reference* CCL HE 12 18
- medium: *reference* CCL HE 12 22
- long: *reference* CCL HE 12 30

Customised abutments

- Implant connections made on 11 axes CNN machines with a precision of 5 µm to guarantee the precision of implant/abutment assembly and sealing.
- Enables perfect defining of the cervical contour.
- Enables treatment of cases that cannot be treated using standard abutments.
- Designed to be compatible with the future tooth and sleeve.



Titanium



Zirconium or
IPS e.max® on
Esthetibase

Simple or anatomical bars and bridges on abutments or direct implants


- High precision.
- Finishing quality.
- Complies with the Sheffield test.



Brand	For suprastructure	Customised abutments		Bars	
		on Esthetibase interfaces	titanium	on implants	on abutments
AB DENTAL®	implant I® platform Ø 3.75 - L compatible	x	x	-	-
	Axiom® compatible	x	x	-	-
ANTHOGYR®	Ossfit® Ø 4.8 compatible	x	x	x (titanium)	x
	Ossfit® Ø 6.5 compatible				
	Anthofit® HE Ø 4.1 compatible	x	x	x	x
ASTRA®	Anthofit® HE Ø 5 compatible				
	3.5 et 4 connection Ocean® compatible	x	x	-	-
BIOHORIZONS®	4.5 et 5 connection Lilas® compatible				
	connection yellow Ø 3.5 - L compatible				
	connection green Ø 4.5 - L compatible	x	x	-	-
	connection blue Ø 5.7 compatible				

Brand	For suprastructure	Customised abutments		Bars	
		on Esthetibase interfaces	titanium	on implants	on abutments
BIOMET 3j®	External Hex. Ø 4.1 mm compatible				
	External Hex. Ø 5 mm compatible				
	External Hex. Ø 6 mm compatible (switching emergence Ø 5 if low angulation)				
	Certain® Ø 3.4 Purple compatible	X	X	X	-
	Certain® Ø 4.1 Blue compatible				
	Certain® Ø 5 Yellow compatible				
	Certain® Ø 6 Green compatible				
	Low Profil - Ø 4.8 compatible	-	-	-	X
BIOTECH®	Kontakt® (Ø 3.6 – 4.2 – 4.8 - 5.4 mm) compatible	X	X	-	X
EASY IMPLANT®	connection OCEAN® compatible	X	X	-	-
	connection LILAS® compatible				
GLOBAL D® (TEKKA® & SERF®)	Tekka® Inkone® compatible	X	X	-	-
	Tekka® Hexa Color® green compatible				
	Tekka® Hexa Color® orange compatible				
	Tekka® Hexa Color® yellow compatible				
	Tekka® Hexa Color® purple compatible	X	X	X	-
	Hexa Color® blue compatible				
	Hexa Color® grey compatible				
	Serf® EVL® Ø 3.3 yellow compatible				
	Serf® EVL® Ø 4 blue compatible				
	Serf® EVL® Ø 5 pink compatible	-	X	X	-
IMPLANT DIRECT®	Legacy® green Ø L compatible				
	Legacy® purple Ø L compatible				
	Legacy® yellow compatible	X	X	-	-
	Interactive® purple Ø S compatible				
	Interactive® yellow Ø L compatible				
MIS®	Seven® M4 SP purple compatible	X	X	-	-
	Seven® M4 WP green compatible				
	Multi Unit Standard Ø 4.8 compatible	-	-	-	X
NOBEL BIOCARE®	Active® Ø 3.5 NP - Replace® CC compatible	X	X	-	-
	Active® Ø 4.3 and 5 RP - Replace® CC compatible				
	Replace® NP red compatible				
	Replace® RP yellow compatible				
	Replace® WP blue compatible	X	X	X	-
	Brånemark® NP compatible				
	Brånemark® RP compatible				
	Brånemark® WP compatible	-	X	X	-
	Multi unit RP compatible	-	-	-	X
STRAUMANN®	Tissue level® RN compatible	X	X	x (titanium)	-
	Tissue level® WN compatible				
	Bone Level® NC Ø 3.3 mm compatible	X	X	-	-
	Bone Level® RC Ø 4.1 and 4.8 mm compatible				
ZIMMER®	Tapered Screw-Vent® Ø 3.5 mm compatible				
	Tapered Screw-Vent® Ø 4.5 mm compatible	X	X	-	-
	Tapered Screw-Vent® Ø 5.7 mm compatible				

x : Available / - : Not available



Part 5

Services

- 188 **Welcome pack**
- 188 **Ordering and organisation**
 - 188 Organisational documents
 - 188 Documents to help with your order
- 189 **Patient information**
 - 189 Patient website
 - 189 Waiting room video
 - 189 Patient card
 - 189 Patient file
 - 189 Patient leaflet
 - 189 Waiting room poster

SERVICES

Welcome pack



Are you a new client?
We will send you a free welcoming pack with your surgical kit, comprising:

- all communication tools to help you communicate with your patients (see list of tools page opposite);
- all support material required to ensure the maintenance of your implants in complete confidence (guide for managing complications, guarantee form etc.);
- all the useful documentation for fitting implants and prosthetic restorations (user guide, slides, monitoring sheet for managing the wear and tear of drills etc.).

Ordering and organisation

Organisational documents



All the tools necessary to streamline surgical procedures:

- Medical questionnaire;
- Informed consent;
- Quote type;
- Statistical summary;
- Model of Xray prescription;
- ...

Documents to help with your order



These documents will help your referrals to choose and order the parts necessary to create the desired prosthetic component. These documents are available for our entire implant range.

Patient information

Patient website



A dedicated website is available to patients to inform them about and the advantages of their implant treatments: www.my-dental-implant.com



Waiting room video



An entertaining video to inform your patients about their implant treatment in your waiting room.

Patient card



A small, pocket-sized leaflet to provide your patients with the exact details of the implants fitted into their mouth for better follow-up.
Pack of 25 copies.

Patient file



Use the patient file for efficient traceability of your implant treatment plans and surgical procedures. The patient file is also a working guide so as not to forget any of the important elements that determine the success of the treatment (pre-implant study, information regarding fitting, procedure steps etc.). The patient file will enable you to ensure follow-up of implant fitting that is as good as that for fitting the prosthesis.
Pack of 25 copies.

Patient leaflet



Leaflet providing the patient with a simple explanation of the implant treatment and its advantages over traditional treatment methods.
Pack of 50 copies.

Waiting room poster



Poster to display in your surgery or waiting room (dimensions: 297 mm x 420 mm).





euroteknika: 726 rue du Général De Gaulle - 74700 Sallanches - France
T: +33 (0)4 50 91 49 20 - sales@etk.dental - www.etk.dental

euroteknika implants are medical devices of Class I b (European Directive 93/42/CEE) comply with the standards of conformity and CE0450 marking center.
Read carefully the instructions for use and user manual.
catalogue_3018_G0_0_1801_web