



IBONE® GUIDED SURGERY TO SECURE YOUR PRACTICE


LYRAETK
DIGITAL DENTAL PROTOCOLS
Simple . Safe . Accessible

100% CUSTOMIZED GUIDED SUPPORT FROM PLANNING TO IMPLANT INSERTION

The kit allows you to insert implants from diameters 3.8 to 6.2 and lengths from 6 mm to 14 mm.

All clinical indications can be performed with this kit:

- From single-unit case, partial to total edentulism
- From immediate functional loading on the framework to deferred loading
- From a conventional prefabricated prosthetic system to a custom prosthetic system through CAD/CAM




CLEAR, FLUID PROTOCOL

The kit includes guided protocol tables and a signposted path that tells you which drill series to use.

4.3	Ø4.8 Ø WP	Ø5.5 Ø WP
D4	1-4+10	1-5
D2 - D3	1-7	1-8
D1	1-7+10	1-9*


A GUIDED DRILL RIGHT FROM THE START

With self-guided drills, including the shortest versions, the complete surgical kit ensures simple and effective guidance.



OPTIMIZED SLEEVES

The complete guided surgery sleeve now features lateral flats for use on recessed teeth. The flats allow a 1 mm reduction in lateral diameter. Milling on the upper part of the sleeve gives an indication of its correct orientation: vestibular side.



A PLACE FOR YOUR IMPLANTS

The kit includes arched locations to place the tubes with the implants in the correct order, in line with the planning to ensure a smooth protocol from bone marking to implant insertion.




DENTAL APP PLATFORM, SIMPLY ORDER AND MANAGE YOUR GUIDES.

Choose the level of support you require.

LYRA ETK PLANS AND PRINTS

1 I submit my impression and my cone beam on the Dental APP.




2 LYRA ETK takes care of implant planning



3 I validate the planning and LYRA ETK prints the guide.

I PLAN, LYRA ETK PRINTS

1 I do the implant planning myself.



2 I upload my guide file (.stl) to the Dental APP.



3 LYRA ETK prints the guide.



4 I perform the surgery.

LYRA ETK PROTOCOLS



- Based on **the most advanced knowledge** of the anatomy and behaviour of peri-implant tissues.
- Less invasive, **more conservative restorations**.
- **Codified protocols**, predictable and reproducible performance.
- **Save time in the dental chair** by optimizing procedures.
- Training in procedures **as close as possible to your own practice** (tutoring) and at your own pace (e-learning).
- **Compatible** with most intra-oral scanners, implant systems, implant planning and design software.

COMMITMENT TO EXCELLENCE

LYRA ETK is committed to excellence and continuous improvement in order to offer the best guarantees.



Implants with lifetime guarantee

Prosthetic parts guaranteed for 10 years

LYRA ETK
DIGITAL DENTAL PROTOCOLS

Simple . Safe . Accessible

LYRA ETK - 726 rue du Général de Gaulle - 74700 Sallanches - France

Service clients France : T. +33 (0)1 76 42 00 60 | International Customer Service: T. +33 (0)1 76 36 14 04
sales@lyraetk.com | www.lyraetk.com